



# Bed Head Unit

## Jiangsu SAIKANG Medical



## Bed Head Unit



Jiangsu SAIKANG Medical Equipment Co.,Ltd.

 : No.35 Lehong Road, Modern Agriculture Demonstration Park, Zhangjiagang City,  
Jiangsu Province, China.

 : 0512-58525850

 : [www.chinasaikang.com](http://www.chinasaikang.com)

 : 0512-58520713

 : [sales@chinasaikang.com](mailto:sales@chinasaikang.com)




## Product manual

Jiangsu SAIKANG Medical Equipment Co.,Ltd.

## MODULAR STRUCTURE



- Adjustable length
- Expandable and upgradeable
- Hygienic surfaces, the surface adopt electrostatic powder coating with static dust,
- Flexible positioning of gas terminal units and electronic connections
- Variety of integrated light options to choose from
- Many model for your choose
- Front panels can be opened, easy to maintenance work
- Security, inside of panel, the unit is divided into several passages by an aluminum barrier.
- extruded aluminum profiles. The length can be freely configured up to a maximum of 6 m.
- For Intensive Care Units, Recovery Rooms and Wards

**Model 1: BH-MEX005**



## Bed Head Unit

**Model 2: BH-CTR001**



### Technical:

**Material:** aluminum alloy

**Dimension:** 183(W) mm x66(D) mm, length defined by customer

Middle panel width is 100mm

**Standard color:** both side is white, middle is gray



**Light:** top and down can install light, for your choose

Medical rail is optional



## Bed Head Unit

**Model 3: BH-GTR001**



### Technical:

**Material:** aluminum alloy

**Dimension:** 195(W) mm x81(D) mm, length defined by customer



Middle panel width is 100mm

**Standard color:** white color

**Light:** top and down can install light, for your choose

## Bed Head Unit

**Model 4: BH-GTE001**



### Technical:

**Material:** aluminum alloy

**Dimension:** 248(W) mm x81(D) mm, length defined by customer

Middle panel width is 100mm




**Standard color:** white color

**Light:** top and down can install light, for your choose

Medical rail is optional

## Bed Head Unit

**Model 5: BH-GTR002**



### Technical:

**Material:** aluminum alloy



**Dimension:** 246(W) mm x89(D) mm, length defined by customer

Middle panel width is 90mm

**Standard color:** white color, yellow color is optional No light

## Bed Head Unit

**Model 6: BH-GTE003**



### Technical:

**Material:** aluminum alloy

**Dimension:** 301(W) mm x81(D) mm, length defined by customer

Middle panel width is 100mm

**Standard color:** white color



**Light:** top and down can install light, for your choose

Medical rail is optional



## Bed Head Unit

**Model 7: BH-MEX003**



### Technical:

**Material:** aluminum alloy

**Dimension:** 244(W) mm x 86(D) mm, length defined by customer

Middle panel width is 100mm


**Standard color:** white color, Middle is gray or white is optional

**Light:** down can install light Medical rail is optional

Monitor Bracket is optional

## Bed Head Unit

**Model 8: BH-MEX001**



### Technical:

**Material:** aluminum alloy



**Dimension:** 176.5(W) mm x 26.3(D) mm, length defined by customer

**Standard color:** white color, yellow is optional

**Light:** down can install light

## Bed Head Unit

**Model 9: BH-MEX002**



### Technical:

**Material:** aluminum alloy

**Dimension:** 211.5(W) mm x 84.6(D) mm, length defined by customer

Middle panel width is 100mm

**Standard color:** both side is white, middle is gray

**Light:** top and down can install light, for your choose

Medical rail is optional

## Bed Head Unit

**Model 10: BH-AWB001**



- Technical:**
- 1、 Push-pull murals can meet the health care functions, and achieve the overall room design at the same time.
  - 2、 The health care equipment at interior, covert design, effectively block dust and man-made pollution.
  - 3、 External mural design, with the overall design of furniture to create comfortable atmosphere for recovery.
  - 4、 Separation of gas and electricity, safe and reliable.
  - 5、 Profuse medical products and accessories, clients can customize