

SAIKANG
MEDICAL

Race against time, healthier your life!



Treatment Table ■■

Addr. No.35 Lehong Road,Zhangjiagang City,Jiangsu Province,China.

Tel. 0086-18021231901

Fax. 0086-512-58525850

Web. <http://www.jssaikang.com>

<http://www.chinasaikang.com>

E-mail sales@chinasaikang.com



JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

Descriptions:

This electric treatment table has multi functions of massage, beauty, tattoo.
PU surface washable mattress with optional upholstery(6cm polyurethane sponge)
Standard hand controller,easy and convenient to control



Face hole with gas spring




4 inch casters



1.5mm thickness powder coated steel structure



6cm PU mattress




Footbar system



Dongrun brand motor

Technical parameter	
Model	X23
Total dimension	1970x690x(530-970)mm
Head section up/down	10°/85°
Leg section up/down	80°/0°
Waist section up/down	30°/0°
Height adjustment	530-970mm
Hand controller	1 pc
Mattress	6 cm PU mattress
Motor brand	Dongrun CE approved motor
Steel	1.5mm powder coated steel
Standard accessories	
4 inch casters	4 pcs
Motor	1 pc
Pillow	1 pc
Gas spring	2 pcs
Hand controller	1 pc
Optional accessories	
Paper roll holder	1 pc
Footbar system	1 pc
Packing	
Packing dimension	200*70*55cm
G.W.	93kgs



Descriptions:

PU surface washable mattress with optional upholstery(6cm polyurethane sponge)
 High quality Virgin ABS material
 Standard hand controller,easy and convenient to control
 The robot welds without scar,standardize laser cutting



Face hole with Rack



6cm PU mattress



2 inch casters



1.2mm thickness steel




Dongrun brand motor




Hand controller

Technical parameter	
Model	X25
Total dimension	1950x680x(600-880)mm
Head section up/down	40°/0°
Leg section up/down	60°/0°
Waist section up/down	20°/0°
Height adjustment	600-880mm
Hand controller	1 pc
Mattress	6 cm PU mattress
Motor brand	Dongrun CE approved motor
Steel	1.2mm powder coated steel
Standard accessories	
2 inch casters	2 pcs
Motor	1 pc
Pillow	1 pc
Rack	2 pcs
Hand controller	1 pc
Optional accessories	
Paper roll holder	1 pc
Footbar system	1 pc
Packing	
Packing dimension	187*71*65cm
G.W.	70kgs



Descriptions:

PU surface washable mattress with optional upholstery(6cm polyurethane sponge)
 High quality Virgin ABS material
 Standard hand controller,easy and convenient to control
 The robot welds without scar,standardize laser cutting



Face hole with Rack



2 inch casters



6cm PU mattress



Hand controller




1.2mm thickness steel



Dongrun brand motor

Technical parameter	
Model	X26
Total dimension	1810x680x(570-860)mm
Head section up/down	50°/0°
Height adjustment	570-860mm
Hand controller	1 pc
Mattress	6 cm PU mattress
Motor brand	Dongrun CE approved motor
Steel	1.2mm powder coated steel
Standard accessories	
2 inch casters	2 pcs
Motor	1 pc
Pillow	1 pc
Rack	2 pcs
Hand controller	1 pc
Optional accessories	
Paper roll holder	1 pc
Footbar system	1 pc
Packing	
Packing dimension	185*69*61cm
G.W.	64kgs



Descriptions:

PU surface washable mattress with optional upholstery(6cm polyurethane sponge)
 High quality Virgin ABS material
 Standard hand controller,easy and convenient to control
 The robot welds without scar,standardize laser cutting




Face hole with gas spring



6cm PU mattress



4 inch casters



Footbar system




1.5mm thickness steel



Dongrun brand motor

Technical parameter	
Model	X27
Total dimension	1970x690x(530-970)mm
Head section up/down	75°/0°
Leg section up/down	80°/0°
Waist section up/down	30°/0°
Height adjustment	530-970mm
Hand controller	1 pc
Mattress	6 cm PU mattress
Motor brand	Dongrun CE approved motor
Steel	1.5mm powder coated steel
Standard accessories	
4 inch casters	4 pcs
Motor	1 pc
Pillow	1 pc
Gas spring	2 pcs
Hand controller	1 pc
Optional accessories	
Paper roll holder	1 pc
Footbar system	1 pc
Packing	
Packing dimension	200*70*55cm
G.W.	93kgs



Descriptions:

PU surface mattress with optional upholstery(6cm polyurethane sponge,10cm manual sponge)
 High quality Virgin ABS material
 Hand controller controls the height adjustment,rack controls the back section adjustment
 The robot welds without scar,standardize laser cutting




6cm PU mattress



2 inch casters



Hand controller



Dongrun hydraulic pump




Face hole with rack



1.2mm thickness steel

Technical parameter	
Model	X22
Total dimension	1970x690x(500-960)mm
Head section up/down	10°/85°
Height adjustment	500-960mm
Foot controller	1 pc
Mattress	6 cm PU mattress
Pump brand	Dongrun CE approved pump
Steel	1.2mm powder coated steel
Standard accessories	
2 inch casters	4 pcs
Hydraulic pump	1 pc
Pillow	1 pc
Rack	2 pcs
Hand controller	1 pc
Optional accessories	
Paper roll holder	1 pc
Packing	
Packing dimension	200*70*55cm
G.W.	73kgs



Descriptions:

PU surface mattress with optional upholstery(6cm polyurethane sponge,10cm manual sponge)
 1.2mm thickness powder coated steel structure
 High quality Virgin ABS material
 Hand controller controls the height adjustment,rack controls the back section adjustment



Face hole with rack



6cm PU mattress



2 inch casters



Hand controller



1.2mm thickness steel



Dongrun hydraulic pump



Technical parameter	
Model	X24
Total dimension	1810x680x(570-860)mm
Head section up/down	50°/0°
Height adjustment	570-860mm
Foot controller	1 pc
Mattress	6 cm PU mattress
Motor brand	Dongrun CE approved motor
Steel	1.2mm powder coated steel
Standard accessories	
2 inch casters	2 pcs
Hydraulic pump	1 pc
Pillow	1 pc
Rack	2 pcs
Hand controller	1 pc
Optional accessories	
Paper roll holder	1 pc
Packing	
Packing dimension	200*70*55cm
G.W.	60kgs

