



JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

Addr. No.35 Lehong Road,Zhangjiagang City,Jiangsu Province,China.
Tel. 0086-18021231901
Fax. 0086-512-58525850
Web. <http://www.jssaikang.com>
<http://www.chinasaikang.com>
E-mail sales@chinasaikang.com



COMPUTER CART



ABOUT US

Jiangsu Saikang Medical Equipment Co.,Ltd.has been established since 2002. It's located in Zhangjiagang City, 2.5 hours' driving from Shanghai city, which has one of the biggest ports in China.

Saikang is a comprehensive medical solution provider which combines with manufacture, sales, service, and supply of Hospital furniture, Emergency and Surgical products. We do our best to create a brilliant future with our customers.

So far, over 90 countries and regions have enjoyed our products and service. As a responsible enterprise, Saikang always follow the corporate philosophy of "Integrity is the base, Customer is the first", and devote ourselves to the mission of "Race against time, healthier your life".

Intelligent Mobile Medication Cart Solutions



New All-in-One Computer Cart	01
New Laptop Cart	02
All-in-One Computer Cart	03
Simple All-in-One Computer Cart A	04
Simple All-in-One Computer Cart B	05
Simple All-in-One Computer Cart C	06
Laptop Cart	07
Simple Laptop Cart	08
Simple Workstation Cart	09

Workstation Cart	10
HD Video Workstation Cart	11
Dossier All-in-One PC Cart	12
Nursing Cart A	13
Nurse Treatment Cart	14
Telemedicine System Cart	15
Bedside Charging Cart	16
Medication Cart A	17
Endoscopic System Cart B	18

ECG Laptop Cart	19
Multi-Functional Small Equipment Cart	20
Small Medical Equipment Cart	21
Elevated Small Medical Equipment Cart	22
Patient Monitor Cart A	23
Patient Monitor Cart B	24





Model: SKR-AB00

YouTube Link: <https://youtu.be/LS7wL34GMww>


- **Optional:** portable plug-in high-performance lithium battery power supply system (with LED keyboard lights), provide uninterruptible power for equipment or computer, safe, reliable and convenient

- **Casters:** dual wheel casters for medical use, mute and smooth, the front two casters equipped with double pedal and total lock, 360 dwegrees brake control, rear followed by two free casters

- **Base:** ABS plastic shell, stability, easy to clean and maintain, highstrength metal structure base frame, strong load capacity and stability, the base integrated with foot pedal to control height adjustment

New All-in-One Computer Cart


- Height adjustment range 250±10mm, loading capacity≥12Kg
- Cart body height: 1145mm—1395mm
- Height of Tabletop: 780mm-1030mm, Two optional control modes, by the foot pedal or the handle control to operate up and down
- Durable ABS plastic tabletop, size: 490mm(L) x 460mm(W), the middle of the table can be opened for placing mini pc
- Optional tabletop material: the top of the table can be embedded in stainless steel panel, wearresistant and easy to clean, strong metallic
- Monitor bracket size: 386mm(L) x 80mm(W) x 64.5mm(H), built-in constant lift structure, light and convenient, display rack hi-low range: 130mm, rack can be tilted 20°, proposed size of display ≤ 22 inch or all-in-one machine, recommended load ≤ 10Kg
- Bi-directional retractable mouse tray, both working by left and right hand, Keyboard tray size: 465mm(L) x 234mm(W)
- The bottom of the tabletop reserved position for built-in battery and inverter,



- Height adjustment range 250±10mm, loading capacity≥12Kg
- Height of Tabletop: 780mm-1030mm, Two optional control modes, by the foot pedal or the handle control to operate up and down
- Durable ABS plastic tabletop, size: 490mm(L) x 460mm(W), the middle of the table can be opened for placing the laptop.Optional tabletop material: the top of the table can be embedded in stainless steel panel, wear-resistant and easy to clean, strong metallic
- Hidden laptop compartment, size: 417mm(L) x 263mm(W) x 50mm(H), can be placed laptop of 17-inch and below, laptop compartment with lock security
- Bi-directional retractable mouse tray, both working by left and right hand, Keyboard tray size: 465mm(L) x 234mm(W)

- **Optional:** Portable plug-in high-performance lithium battery power supply system (with LED keyboard lights), provid e uninterruptible power for equipment or computer, safe, reliable and convenient
- **Casters:** dual wheel casters for medical use, mute and smooth, the front two casters equipped with double pedal and total lock, 360 degrees brake control, rear followed by two free casters
- **Base:** ABS plastic shell, stability, easy to clean and maintain, high-strength metal structure base frame, strong load capacity and stability, the base integrated with foot pedal to control height adjustment

Model: SKR-BB00



New Laptop Cart




All-in-One Computer Cart




- Base: PA6 + ABS cover
- Battery and UPS system inside the base; White LEDs provide light over keyboard tray
- Front-locking silent casters

- Height adjustment range for LCD and worksurface: 400mm±20mm, Lifting load ≤ 15 kg
- Worksurface height: 890mm-1290mm
- Large, smooth ABS worksurface dimensions: 533mm(L) x 455mm (W), Proposal Table load ≤ 5 kg
- Drawer inner dimensions: 379mm(L) x 298.5mm(W) x 84mm(H)
- Keyboard tray dimensions: 440mm(L) x 256mm(W)
- Mouse tray dimension: 208mm(L) x 200mm(W)
- Storage Bin with handle
- LCD column dimensions: 386mm(L) x 80mm(W) x 64.5mm(H)
- Height adjustable with the help of internal lift construction.
- LCD hanger: height adjustment for LCD pivot is 13cm; 20° tilting for LCD; Applied for ≤ 22 inch LCD display.
- Proposal load ≤ 10 kg



Model: SKR-A00

Simple All-in-One Computer Cart A



Model: SKR-A01


- Concise and compact for 22" All-in-One computer, Proposal load ≤ 12 kg
- Keyboard/mouse tray height: 776mm
- Keyboard tray dimensions: 440mm(L) x 256mm(W)
- Mouse tray dimension: 208mm(L) x 200mm(W)
- Back storage bin with handle, handle height: 790.6mm
- Printer tray dimensions: 310mm(L) x 199mm(W), height 297.6mm
- Cart dimensions: 648.3mm(L) x 489.5mm(W) x 1260.5mm(H), Built in power line
- Built in power line
- Silent medical casters



Simple All-in-One Computer Cart B



YouTube Link: <https://youtu.be/QMPlgwysK40>



Model: SKR-B00

- Concise and compact,height adjustment range for worksurface:400±20mm,lifting load ≤ 15Kg
- LCD column dimesion: 386mm(L) X 80mm(W) X 64.5mm(H),height adjustable with the help of internal lift construction.Applied for 22 inch LCD display, proposal load ≤ 12Kg
- Large,smooth ABS worksurface dimension:533mm(L) X 455mm(W), worksurface height:750mm-1100mm, Can be customized built-in mini computer host
- Optional rotary keyboard tray, keyboard tray dimension:440mm(L) X 256mm(W)
- Mouse tray dimension:208mm(L) X 200mm(W)
- Storage Bin with handle
- Optional Built in power line
- Silent medical casters

Model: SKR-B00P



PC + Printer


Mini PC + Printer

Simple All-in-One Computer Cart C



Model:SKR-C00

- Concise and compact,gas spring lifting column, pedal plate lift and adjust,height adjustment range for worksurface:230±10mm,lifting load ≤ 18Kg
- LCD column dimesion: 386mm(L) X 80mm(W) X 64.5mm(H),height adjustable with the help of internal lift construction.Applied for 22 inch LCD display, proposal load ≤ 10Kg
- Metal worksurface dimension:470mm(L) X 352mm(W) or 470mm(L) X 460mm(W),worksurface height:898mm-1128mm or 933mm-1163mm
- Retractable metal keyboard tray, keyboard tray dimension:470mm(L) X 206mm(W)
- The rotary type mouse tray , mouse tray dimension:189mm(L) X 195mm(W)
- Under the table can be built into the 2 sets of plug-in lithium power supply module, optional inverter or direct output, fully meet the needs of mobile information
- Base: SPCC (steel plate) + ABS cover
- Silent medical casters



Model: SKR-C00H



lithium power supply module



Mini PC module

Laptop Cart




Model:SKR-B10

Model: SKR-B12




- Height adjustment range for worksurface: 400±20mm, Lifting load ≤ 15kg, Worksurface height 890mm-1290mm
- Large, smooth ABS worksurface dimensions:533mm(L) x 455mm (W)
- Laptop cabinet dimensions: 430mm(L) x 346mm(W) x 56.5mm(H)
- With external 2-Outlet power strip and main power cable at the cart base.
- Applied for laptop(≤ 15")
- Internal dimensions for drawer: 379mm(L) x 298.5mm(W) x 84mm(H)
- Keyboard tray dimensions: 440mm(L) x 256mm(W)
- Mouse tray dimension: 208mm(L) x 200mm(W)
- Base: PA6 + ABS cover
- Silent medical casters



Simple Laptop Cart



- Height adjustment range for worksurface:400 ±20mm,worksurface height:750mm-1150mm,lifting load ≤ 15Kg
- Large,smooth ABS worksurface dimension:533mm(L)×455mm(W)
- Laptop cabinet dimension:430mm(L)×346mm(W)×56.5mm(H)
- Applied for ≤ 15 inch laptop
- Optional rotary keyboard tray, keyboard tray dimension:440mm(L)×256mm(W)
- Mouse tray dimension:208mm(L)×200mm(W)
- Base:PA6 + ABS cover




Model: SKR-IB00

YouTube Link: <https://youtu.be/EZOugngltUo>

Simple Workstation Cart



Model: SKR-J00




- Simple and practical, the cart with PC tray and printer tray for all types of workstations PC and printers
- The integrated design of the mouse bracket and the display bracket, applied for ≤ 24 inch LCD display, proposal load ≤ 15 Kg
- Keyboard tray dimensions: 470mm(L) x 206mm(W)
- The rotary type mouse tray, dimensions: 189mm(L) x 195mm(W), with mouse rack
- Behind of the LCD fix frame is handle, easy to move the trolley
- Base: SPCC (steel plate)
- Cart dimensions(not include equipment): 706.5mm(L) x 572mm(W) x 1173(H)
- Silent medical casters

Workstation Cart



Model: SKR-G20




- Big capacity box for lots of workstation device
- Large, smooth ABS worksurface dimensions: 533mm(L) x 455mm (W), Proposal Table load ≤ 10 kg
- LCD column dimesion: 386mm(L)x80mm(W)x64.5mm(H), height adjustable with the help of internal lift construction. Applied for 22 inch LCD display, proposal load ≤ 8 Kg
- Worksurface height: 960 mm
- Built in power line
- Keyboard and mouse tray, Back storage bin with handle
- Cart dimensions: 535mm(L) x 650mm(W) x 1350(H)
- Silent medical casters

HD Video Workstation Cart



Model:SKR-G21

- **Optional:** Battery and UPS System
- Silent medical casters
- Built-in computer host specification:
CPU: Intel E3 1231 V3
Mother board: Asus B85 PLUS
Memory: 8G/1600
Graphics card: DDR5 2G
HD: WD 1000G/SATA
Video Capture Card




- Height adjustment range for LCD and worksurface: 400mm
- Worksurface height: 890mm-1290mm, Lifting load ≤ 15 kg
- ABS worksurface dimensions: 533mm(L) x 455mm(W)
- Keyboard tray dimensions: 440mm(L) x 256mm(W)
- Mouse tray dimension: 208mm(L) x 200mm(W)
- LCD column dimensions: 386mm(L) x 80mm(W) x 64.5mm(H)
- Height adjustable with the help of internal lift construction. Proposal load ≤ 10 kg
- LCD hanger:height adjustment for LCD pivot is 13cm; 20° tilting for LCD for ergonomic comfort
- Base: PA6 + ABS cover

Dossier All-in-One PC Cart



Model:SKR-G00


- Battery and UPS system inside the base;White LEDs provide light over keyboard tray;used for ≤ 22 inch All-in-One PC . Proposal load ≤ 10 kg
- 20 tiers inside dossier,tier height 20mm,width 240mm
- Height from countertop to floor: 1048.5mm, Proposal Table load ≤ 15 kg
- Large, smooth ABS worksurface dimensions:533mm(L) x 455mm (W)
- Ergonomic back-tilt keyboard tray with left or right mouse tray,height:828.4mm
- Keyboard tray dimensions: 440mm(L) x 256mm(W)
- Mouse tray dimension: 208mm(L) x 200mm(W)
- Back storage Bin with handle,handle height: 830mm
- Cart dimensions: 659.1mm(L) x 640.9mm(W) x 1436.6mm (H)
- Front-locking silent casters



Nursing Cart A



Model:SKR-NC1




Model: SKR-NC2

- Cart Body made of aluminum alloy with steel box frame structure, light and strong, also corrosion resistance
- Six spacious drawers,5 drawer with different size of lattices for different medicine. Drawers by interlocking system.
- Worksurface height:977mm,10mm thickness tempered glass worksurface, easy-to-clean,dimensions: 500mm(L) x 500mm(W), Proposal Table load ≦ 5kg- Left side with sterilize liquid bottle box and document basket
- Right side with 2 sharps boxes and 2 garbage bins
- Optional Laptop or all-in-one computer rocker arm and scanner gun hanger
- 768mm(L)×585mm(W)×1377mm(H) (include LCD arm height)
- Silent medical casters
- Option :Color Pink or Blue. Retractable stainless steel infusion stand , lithium power system

Nurse Treatment Cart

Model:SKR-DB5

- Height adjustment range for worksurface: 400±20mm, Lifting load≤15kg
- Worksurface height:750-1170mm (±20mm)
- Worksurface dimensions: 610mm(L) x 470mm(W)
- There is a small laptop(≤14") cabinet which is lockable with 2 keys
- 2 big storage trays with different depth.
- Dimensions: 469.5mm(L) x 306mm(W)
- Fence with hangers(include 3 small pockets,a scanner gun holder, a gel cup)
- Base: PA6 + ABS cover
- Silent medical casters



Model: SKR-J205


Model: SKR-D05

360 °wheels



Telemedicine System Cart

- Height adjustment range for LCD and worksurface: 400±20mm
- Total height is 1534mm-1934mm, Lifting load ≤ 17kg
- Applied for 2 ≤ 21" LCD, Proposal evaluation load ≤ 12kg
- HD 360 degree rotating webcam with remote controller
- Worksurface height: 890-1290mm
- ABS worksurface dimensions: 533mm(L) x 455mm(W)
- Keyboard tray dimensions: 440mm(L) x 256mm(W)
- Mouse tray dimension: 208mm(L) x 200mm(W)
- Silent medical casters
- Back storage bin with handle,handle height: 680mm-1080mm
- Built-in computer host specification(for reference) :
CPU: Intel E3 1231 V3




- **Optional:** Battery and UPS System
- **Mother board:** Asus B85 PLUS
- **Memory:** 8G/1600
- **Graphics card:** DDR5 2G
- **HD:** WD 1000G/SATA
- Video Capture Card

Model:SKR-F00



Bedside Charging Cart

- Convenient for bedside charging in sickroom
- Worksurface height: 900mm
- Arm for all-in-one computer or laptop
- Drawer space for POS device, keyboard, mouse
- Counting machine tray
- Optional document basket, grid basket
- Cart dimensions: 1000mm(L) x 510mm(W) x 980mm(H) (not include LCD arm height)
- Battery and UPS system
- Silent medical casters




Model:SKR-BCC1



Medication Cart A



Model:SKR-CB1



- Six spacious drawers,5 drawer with different size of lattices for different medicine.
- Drawers by interlocking system.
- Worksurface height: 990mm
- 10mm thickness tempered glass worksurface,easy-to-clean dimensions: 498mm(L) x 500mm(W)
- Proposal Table load \leq 5kg
- Foldable worksurface dimensions: 470mm(L) x 290mm(W)
- Laptop rocker arm and scanner gun hanger
- 2 Dustbins
- Cart dimensions: 700mm(L1)/950mm(L2) x 600mm(W) x 1200mm(H) (include rocker arm height)
- Silent medical casters

Endoscopic System Cart B




Model:SKR-ESC2




Model: SKR-ND2

- Convenient for endoscopic system device
- Layers and layer height adjustable, Proposal load \leq 20kg/Layer
- Cart dimensions: 664mm(L) x 630mm(W) x 1650(H)
- Worksurface height: 1125mm
- Worksurface dimensions: 527mm(L) x 420mm(W)
- Rocker arm for display,height-adjustable,360°rotation for LCD, Standard Load \leq 9Kg, optional maximum load \leq 21Kg
- Special power wiring slot
- Silent medical casters



ECG Laptop Cart



Model:SKR-EB5


- It is for mini ECG device with wire sections arm.
- Height adjustment range: 400±20mm. Worksurface height: 890mm-1290mm Lifting load \leq 15kg
- Back cover clamshell design, use for adapter, Built in power line
- Large, smooth ABS worksurface dimensions:533mm(L) x 455mm (W)
- Sliding keyboard/mouse tray height: 684mm-1084mm
- Keyboard tray dimensions: 440mm(L) x 256mm(W)
- Mouse tray dimension: 208mm(L) x 200mm(W)
- Storage bin with handle
- Wheel-to-Wheel Base Dimensions: 513mm x 672.8mm
- Silent medical casters

Multi-Functional Small Equipment Cart




Model:SKR-KB0

- The five-star type base, stable and applied for different small equipment
- Pipe/Wire holder arm available
- Proposal Table load \leq 30kg, proposal mounted load \leq 15kg
- Mesh basket for placement of related goods
- handle height: 911mm
- Metal keyboard/mouse tray,ABS mouse rack, the tray height:818.8mm
- The five-star type base dimensions: The diameter of 450mm
- Silent medical casters



Small Medical Equipment Cart




Model:SKR-KB1



- Applied for different small equipment
- Equipment tray dimensions: 310mm(L) x 235mm(W)
- Pipe/Wire holder arm available
- Equipment tray height: 1012mm, Proposal Table load $\leq 30\text{kg}$
- Mesh basket dimensions: 400mm(L)x220mm(W1) / 130mm(W2)x100mm(H)
- Storage bin with handle,handle height: 929.2mm
- Cart dimensions (not include equipment): 565.1mm(L) x 461.8mm(W) x 929.2mm(H)
- Silent medical casters

Elevated Small Medical Equipment Cart




Model:SKR-K00



- The whole design is simple and beautiful, gas spring lifting column, pedal plate lift and adjust
- Can be customized according to customer demand to placed small devices or display equipment
- The tray height: 869mm-1099mm, Proposal Table load $\leq 12\text{kg}$
- Mesh basket dimensions: 380mm(L) x 130mm(W) x 100mm(H)
- Standard configuration: Metal bracket of MINI host and keyboard , display rack and push handle
- Base: SPCC (steel plate)
- Standard Cart dimensions (not include equipment): 470mm(L) x 571mm(W) x 1153mm-1383mm(H)
- Silent medical casters

Patient Monitor Cart A



- Concise and compact for 13" patient monitor , Proposal Table load \leq 15kg
- The new cart integration aluminum alloy column, the top equipment mesa, column mounted attachments, easy to transport and install
- The equipment mesa integrated with the handle , the mesa height: 1151.3mm
- The free PABS tray , dimensions: 388mm(L) x 238mm(W) x 50mm(H)
- ABS tray height: 795.9mm
- Cart dimensions: 502.2mm(L) x 629.2mm(W) x 1151.3mm(H)
- 4PCS of 75mm silent medical casters

Model:SKR-J21




Model:SKR-J10



Patient Monitor Cart B



Model:SKR-J20



- Concise,light and compact for PAD patient monitor , Proposal Table load \leq 5kg
- The 40mm stainless steel tubular column , The top is an integrated design of the handrail and the equipment platform , easy to install and use ,
- The equipment mesa height: 839.5mm (Can be customized according to height or lift demand)
- The basket dimensions: 280mm(L) x 130mm(W) x 100mm(H)
- Five foot metal base, the rotation is flexible, stable and firm
- Cart dimensions: 558.3mm(L) x 531.8mm(W) x 999.2mm(H)
- 4PCS of 75mm silent medical casters

Model:SKR-J22

