



SK-P009

Double Arm Electric Surgical pendant



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

Use: ICU room, Suction ceiling installation.

Material: Main body use aluminium alloy material.
Surface use powder-coated treatment.



SK-P001

ICU Pedant Ceiling-Mounted Rail System(Cantilever)

Configurations:

SK-P001

Gas outlets standard	Oxygen	1 pc
	Vacuum	1 pc
	Medical Air	1 pc
Lighting lamps		2 pcs
Swing-arm lamp		1 pc
Electric sockets		6 pcs
Ground terminal		2 pcs
Net connect		RJ45 1 pc
Trays		3 pcs
Drawer		1 pc


Parameters:

SK-P001

Working Power	Customerization
Bridge length	2500-3200mm
Left and right move distance	300-600mm
Horizontal rotation angle	≤340°
Rotation angle	0-340°
Loading capacity	≤150kg
Gas outlets life	Insert&pull out more than 20000times

Use: ICU room, Suction ceiling installation.

Material: Main body use aluminium alloy material.
Surface use powder-coated treatment.



SK-P003

ICU Pedant Ceiling-Mounted Rail System(Fixed Rotary)



SK-P002

Configurations:

SK-P002/SK-P003

Gas outlets standard	Oxygen	1 pc
	Vacuum	1 pc
	Medical Air	1 pc
Lighting lamps		2 pcs
Swing-arm lamp		1 pc
Electric sockets		6 pcs
Ground terminal		2 pcs
Net connect		RJ45 1 pc
Trays		3 pcs
Drawer		1 pc
Adjustable stainless steel iv pole		
Syringe pump support		1 pc
Storage basket		1 pc


Parameters:

SK-P002/SK-P003

Working Power	Customerization
Bridge length	2500-3200mm
Left and right move distance	300-600mm
Horizontal rotation angle	≤340°
Rotation angle	0-340°
Loading capacity	≤150kg
Gas outlets life	Insert&pull out more than 20000times


Use: ICU room, Suction ceiling installation.

Material: Main body use aluminium alloy material.
Surface use powder-coated treatment.



SK-P004

ICU Pedant Ceiling-Mounted Rail System



SK-P005

Configurations: SK-P004

Gas terminal	Oxygen	1 pc
	Vacuum	1 pc
	Medical Air	1 pc
Instrument Platform		4 layers
Swing-arm lamp		1 pc
Lighting Lamps		2 pcs
Ground terminal		2 pcs
Power outlet		6 pcs
Network Interface		RJ45 1 pc
Drawer		1 pc

Parameters: SK-P004

Working Power	Customerization
Bridge length	1800-2200mm
Arm length	400mm
Horizontal rotation angle	$\leq 340^\circ$
Loading capacity	$\leq 150\text{kg}$
Gas outlets life	Insert&pull out more than 20000times

Use: Operating room, Suction ceiling installation.

Material: Use aluminium alloy material
Surface use powder-coated treatment



SK-P008

Double Arm Electric Surgical pendant



SK-P009

Configurations: SK-P008

Gas connecting	Oxygen	1 pc
	Vacuum	1 pc
	Medical Air	1 pc
Rotating horizontal arm		1 pc
Electrical outlet		6 pcs
Gas braking system		1 set
Ground terminal		1 set
Adjustable stainless steel IV pole		1 pc
Adjusted ss injection pump shelf		1 pc

Parameters: SK-P008

Working Power	Customerization
One pc	
With three layers instrument platform	
Horizontal rotation angle	$\leq 340^{\circ}$
Loading capacity	$\leq 150\text{kg}$
Gas outlets life	Insert&pull out more than 20000times

Use: General use, Suction ceiling installation.

Material: Use aluminium alloy material
Surface use powder-coated treatment



SK-P010

Medical Multi-function Column



SK-P011

Medical Floor Stand Column

Configurations:

SK-P010

Gas outlet standard	Oxygen	2 pcs
	Vacuum	2 pcs
	Medical Air	2 pcs
Electric sockets		12 pcs
Adjustable stainless steel iv pole		1 pc
Layers		3 pcs
Syringe pump support		1 set
Drawer		1 pc
Storage basket		2 pcs

Parameters:

SK-P010

Working Power	Customerization
Horizontal rotation angle	≤340°
Loading capacity	200KG
Round corner anti-collision design	
Gas outlets life	Insert&pull out more than 20000times

Use: Operating room, ceiling type.

Material: Use aluminium alloy material
Surface use powder-coated treatment



SK-P014

Single Arm Surgical Pendant(Adjustable Height)



SK-P013

Single Arm Electric Surgical Pendant



SK-P018

Single Arm Surgical Pendant (Adjustable Height)

Configurations:		SK-P014
Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
Plug		8 pcs
Ground terminal		2 pcs
Net work port RJ45		1 set
Adjustable ss IV pole		1 pc
Instrument platform		2 layers
Syringe pump support		1 set
Drawer		1 pc

Parameters:	SK-P014
Working Power	Customerization
Working platform vertical lift range	600mm
Horizontal arm moving range	≤1000mm
Horizontal rotation angle	≤340°
Gas outlets life	Insert&pull out more than 20000times

Use: Operating room, ceiling type.

Material: Use aluminium alloy material

Surface use powder-coated treatment



SK-P012

Double Arm Surgical Pendant



SK-P016

Double Arm Electric Surgical pendant
(height Adjustable)

Configurations:		SK-P012
Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
Electricity socket		8 pcs
Ground terminal		2 pcs
Instrument platform layers		2 pcs
Net work port RJ45		1 set
Adjustable ss iv pole		1 pc
Syringe pump support		1 set
Drawer		1 pc

Parameters:	SK-P012
Working Power	Customerization
Workbench vertical lift	600mm
Radius of cross arm	1400mm
Horizontal rotation angle	≤340°
Gas outlets life	Insert&pull out more than 20000times

Use: Abdominal Cavity , Ceiling type.

Material: Use aluminium alloy material
Surface use powder-coated treatment



SK-P020

Single Arm Electric Abdominal Cavity Pendant (Adjustable Height)



SK-P021

Double Arm Abdominal Cavity Pendant
(Fixed Height)

Configurations:

SK-P020

Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
	Air	1 pc
	CO2	1 pc
Electrical outlet		8 pcs
Ground terminal		2 pcs
Internet connector RJ45		1 pc
Adjustable stainless steel iv pole		1 pc
Syringe pump support		1 set
Drawer		1 pc
Syringe pump support		1 set
Drawer		1 pc

Parameters:

SK-P020

Working Power	Customerization
Platform vertical lift range	600mm
Horizontal arm moving range	≤1000mm
Horizontal rotation angle	≤340°
With gas braking system	
Loading capacity	1400mm
With two layers instrument platform	
Gas outlet insert and draw life	20000 times

Use: Anesthesia, Ceiling type.

Material: Use aluminium alloy material
Surface use powder-coated treatment



SK-P023

Single Arm Electric Anesthesia Pendant (Electric Height Adjustable)

Configurations:		SK-P023
Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
Electric Sockets		8 pcs
Ground terminal		2 pcs
Net work port RJ45		1 pc
Instrument platform		2 layers
Suspension anesthesia machine		
Adjustable ss iv pole		1 pc
Syringe pump support		1 set

Parameters:	SK-P023
Working Power	Customerization
Workbench vertical lift	600mm
Lifting speed	22mm/s
Radius of the cross arm	≤1000mm
Horizontal rotation angle	≤340°
Gas outle Insert&Draw Frequency	20000 times

Use: Endoscope, Ceiling type.

Material: Use aluminium alloy material
Surface use powder-coated treatment



SK-P024

Rotary Endoscope Pendant (Electric Height Adjustable)

Configurations:		SK-P024
Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
Electrical outlet		8 pcs
With equipotential ground terminal		2 pcs
Internet connector RJ45		1 pc
Adjustable stainless steel iv pole		1 pc
Syringe pump support		1 set
Drawer		1 pc

Parameters:	SK-P024
Working Power	Customerization
Horizontal rotation angle	≤340°
With gas braking system	
Loading capacity	100-200KGS
With three layers instrument platform	

Use: Endoscope, Ceiling type.

Material: Use aluminium alloy material
Surface use powder-coated treatment



SK-P025

Rotary Endoscope Pendant (Fixed Height)

Configurations:		SK-P025
Gas terminal	Oxygen	2 pcs
	Vacuum	2 pcs
Instrument platform		2 layers
Electric socket		8 pcs
Grounding terminal		2 pcs
Adjustable ss iv pole		1 pc
Syringe pump support		1 set
Drawer		1 pc
Anti-collision round corner design		

Parameters:	SK-P025
Working power	AC220V, 50HZ
Workbench vertical range	300mm
Electric power	1 kw
Lifting speed	22 mm/s
Horizontal rotation angle	≤340°
Loading capacity	100kg
Gas outle Insert & Draw Frequency	20000 times

Use: ICU/Infusion room, Suction ceiling installation.



SK-P007

Infusion Rails System

Parameter : (Price/Meter)
 Working Power: Customization
 Loading capacity: $\leq 150\text{kg}$
 Different color outlet prevent connect wrong
 Insert and pull out more than 20000times
 Electrical socket: according to drawbridge length
 Infusion hook: according to drawbridge length



SK-P026

Rotary Surgical Pendant (Electric Height Adjustable)



SK-P027

Rotary Surgical Pendant (Fixed Height)