


Addr.: No.35 Lehong Road, Modern Agriculture Demonstration Park,
 Zhangjiagang City, Jiangsu Province, China.
 Tel.: +86 512 58525850
 Fax: +86 512 58520713
 Mob.: +86 13962479538
 P.C.: 215600
 Web.: www.chinasaikang.com www.jssaikang.com
 E-mail: inquiry@jssaikang.com



About Us

Jiangsu Saikang Medical Equipment Co., Ltd. has been established since 2002. It is located in Zhangjiagang City, along the Economic circle of Yangtse River Delta, enjoying the highly developed economic expansion and convenient transportation access.



Our Workshop



15
years



400
employees



30,000
square meters



Our Team



So far, over 90 countries and regions have enjoyed our products and service. Our success has been mainly achieved by following the corporate philosophy of being firmly committed to "Integrity is the base, Customer is the first."

Advantage of Saikang patient monitor



EM Series Integrated healthcare monitor



Models: EM032 12 for 12"
EM032 10 for 10"
EM032 8 for 8"
Bed side/Healthcare department/Clinical

Basic configuration:
EGC, RESP, HR, SPO2, PR, NIBP, TEMP


Optional configuration:
Touch screen, 12 lead EGC, WIFI&GPRS, 2-IBP, ETCO2, 5AH battery, Central monitoring system, Wall standing & Trolley

Feature

12.1"/10.4"/8.4" high resolution TFT, color screen
Fanless design
ST and arrhythmia analysis
SPO2 low perfusion index up to 0.1%
Three alarm lights well distinguish technical and physiological alarm
Pluggable Li-ion battery, max 5 hours continuous working time Connected with CMS by wire and wireless (3G and WIFI)
Multi-language



Patient Monitoring(ICU/CCU/OR)



Models: EM033 12S for 12"
EM033 10S for 10"
EM033 8S for 8"

ICU/CCU/OR



Advantage of Saikang patient monitor



Feature

- 12.1"/10.4"/8.4" high resolution TFT, color screen
- Semi-Modular and fanless design
- Multi-language
- Various of operation modes: touch screen, knob, mouse, keyboard
- 3/5 lead ECG
- ST and arrhythmia analysis
- STO2 low perfusion index up to 0.1%
- Three alarm lights well distinguish technical and physiological alarm
- Modular printer with 3 waveform output
- Pluggable Li-ion battery, max 5 hours continuous working time
- Connected with CMS by wire and wireless(3G and WIFI)



Basic Configuration:
ECG, RESP, HR, SPO2, PR, NIBP, TEMP

Optional Configuration:
Touch Screen, Recorder, 12 lead ECG, ETO2,
2-IBP, WIFI&GPRS, 5AH battery, Wall standing
&Trolley, Central monitoring system



Central Monitoring System

Single/dual screen monitoring
Monitoring 64 beds
Wire and wireless(3G and WIFI)CMS
Support HL7, HIS



Pad Monitor

First aid, clinic examination, ward round, healthcare center, nursing homes, family healthcare



Standard configuration:

Spo2/PR, ECG/HR, NIBP, Infrared forehead temperature,

Optional configuration:

Blood glucose, CHOL(cholesterol), UA(uric acid)

Feature

Android system, easy operation as iPad
Small size & light weight
10.1" TFT color touch screen
Support Blood glucose, Forehead TEMP, Cholesterol, Uric acid.
Rechargeable lithium battery, 10 hours continuous working time 120 hours data memory

