

Electric Operating Table



Parameters	Specification	Standard accessories	Qty
Length	1960mm	Anaesthetic Screen	1
Width	500mm	Shoulder Rest	2
Height	690-1100mm	Arm Support	2
		Body Support	2
		Knee Crutch	2
		Matress	1



Optional traction 1005 supine position



Optional 1005 lateral lying position



Optional 1005 Tibia extension position



Joint for 1005



Optional 1005 Tibia bone-nailing position



Optional traction 1006



Optional spondyle operating holder



Optional lateral arm support



Anorectum operation position



Optional basin for gynecology



E.N.T. Operation position



Lateral lying position for chest, kidney,cholecyst operation

Electro-hydraulic operating table A106-3

Introduction:

Electro-hydraulic operating table is made specially to meet demand of general surgical operations. It adopts full 304 stainless steel.Tabletop is X-raying available.Auto-lock on remote control, electric base lock and release.Sliding can be 350mm,0 reset function,with battery available.

Technical Data:

Trendelenburg: 30°
Reverse Trendelenburg: 30°
Lateral Tilt(left/right): 18°/18°
Head Rest(up/down): 40°/90°
Footplate(up/down&horizontal positioning): 15°/90° 90°
Back plate(up/down): 90°/40°
Longitudinal Sliding: 350mm
Power: 220v, 50Hz, 1.0kW

Electric Operating Table



Electric Operating Table A106-2

Name: SKF-C Electric Operating Table
Introduction:
Electric operation table is made specially to meet demand of sitting posture operations with a extra low position.It adopts full 304 stainless steel,excellent properties on anti-rust.Tabletop is X-raying available.

Technical Data:
Trendelenburg: ≥22°
Reverse Trendelenburg: ≥15°
Lateral Tilt(left/right): 15°/15°
Head Rest(up/down): 45°/90°
Footplate(up/down&horizontal positioning): 15°/90° 90°
Back Plate(up/down): 75°/20°
Power: 220v, 50Hz, 1.0Kw

Parameters	Specification	Standard accessories	Qty
Length	1960mm	Anaesthetic Screen	1
Width	480mm	Shoulder Rest	2
Height	500-800mm	Arm Support	2
		Body Support	2
		Knee Crutch	2
		Matress	1



Tabletop headward sliding / toward leg sliding



Tabletop elevation up / down



Lateral tilt left / right



Trendelenburg / anti-trendelenburg



Tabletop flexural position



Brake system lock / release



Electric Operating Table




Manual Hydraulic Operating Table A3001-6

Name: SKF-C Manual Hydraulic Operating Table

Introduction:
It is all mechanically operated, hydraulic pressure increase or decrease. It adopts stainless steel as its material with nice appearance and compact structure, made of full 304 stainless steel. Germany made Gas spring to adjust the backrest and leg plate. Tabletop is X-raying Available.

Technical Data:
Length: 2000mm
Width: 500mm
Elevation(up/down): 705mm-1005mm
Trendelenburg/Reverse Trendelenburg: 30°/25°
Lateral Tilt: 20°/20°
Head rest(up/down): 50°/60°
Leg Plate(up/down): 15°/ 90°
Back Plate(up/down): 70°/30°

Parameters	Specification	Standard accessories	Qty
Length	2000mm	Anaesthetic Screen	1
Width	500mm	Shoulder Rest	2
Height	705-1005mm	Arm Support	2
		Body Support	2
		Knee Crutch	2
		Cushion	1




Electric operating table A100-4A

Name: SKF-C Electric operating table

Features:
Table height is adjustable by motor, from 710-1010 mm
Whole table can realize reversed trendelenburg(25°) and trendelenburg (25°)
Whole table can tilt right/left 15°/15°
Head board can be upward or downward , upward 0-40°, downward 0-90°
Back board can be upward or downward , upward 0-90°, downward 0-15°
Kidney bridge elevation range:0-110mm
Horizontal sliding: 350mm
Leg boards can downward 90°, also can be open by 90° to be convenient for medical staff to operate or to be apart according to requirements

Parameters	Specification	Standard accessories	Qty
Length	2060mm	Bed frame	1
Width	500mm	Matress	1
Height	710-1010mm	Anaesthetic Screen	1
		Shoulder Support with cushion	2
		Arm support with cushion	2
		Body holder	2
		Kidney bridge handle	1
		Leg support with cushion	2
		Hand controller	1



Electric Operating Table



Low-position: 750mm



High-position: 50mm



Kidney bridge elevation: 0-80mm



Optional Head Rest



Optional Lateral Arm Support



Optional Waste Basin for Gynecology



Reversed Trendelenburg: 0-22°



Trendelenburg: 0-25°



Anorectum Operation Position



Longitudinal adjustment: 0-295mm



Lateral Tilt Right: 0-20°



Lateral Tilt Left: 0-20°



Optional Orthopedics Traction Frame



Lateral Tilt Right: 0-20°



Lateral Tilt Left: 0-20°

Electric Operating Table A100-4

Name: SKF-C Electric operating table

Introduction:

Electric operating table is made specially to meet demand of general surgical operations. It adopts full 304 stainless steel, excellent properties on anti-rust.

Tabletop is X-raying available, can be sliding for over 300mm, that meet C-arm need.

Technical Data:

- Trendelenburg: $\geq 25^\circ$
- Reverse Trendelenburg: $\geq 25^\circ$
- Lateral Tilt(left/right): $15^\circ/15^\circ$
- Head Rest(up/down): $45^\circ/90^\circ$
- Footplate(up/down&horizontal positioning): 90°
- Back plate(up/down): $75^\circ/20^\circ$
- Kidney bridge: 120mm
- Sliding: 300mm
- Power: 220v, 50Hz, 1.0kW

Parameters	Specification	Standard accessories	Qty
Length	2070mm	Anaesthetic Screen	1
Width	550mm	Shoulder Rest	2
Height	700-1000mm	Arm Support	2
		Body Support	2
		Knee Crutch	2
		Matress	1

Electric Obstetric Bed



Electric Obstetric Bed A99-5

Name: SKF-C Electric Obstetric Bed

Introduction:

Electric operation table for gynaecology and obstetrics. With full affinity profile. Any colors are available for the mattress. Made of full premium 304 stainless steel. Double control(hand remote control and pedal switch)

Technical Data:

Lenth: 1780mm
Width: 600mm
Elevation(up/down): 680mm-880mm
Seat Plate(up/down): 17°/ 6°
Back Plate(up/down): 80°/20°
Power: 220v, 50Hz, 1.0kw

Parameters	Specification	Standard accessories	Qty
Length	1780mm	Anaesthetic Screen	1
Width	600mm	Shoulder Rest	2
Height	680-880mm	Arm Support	2
		Body Support	2
		Knee Crutch	2
		Matress	1

