

24hrs: 0086-18021231901
www.chinasaikang.com

SAIKANG (USA Branch)
Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (East Africa Branch)
Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (China Headquarter)
No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@ssaikang.com

SAIKANG (UK Branch)
Uk Healthward International Limited
Owner(s) Address DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (Mid East Branch)
Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Factory Tour

Website



SAIKANG™

LEADING BRAND OF CHINA
MEDICAL FURNITURE

ABS | PHENOLIC | COLOR PAINTED | STAINLESS | STEEL | SERIES

JIANGSU SAIKANG MEDICAL EQUIPMENT CO., LTD.
www.chinasaikang.com

Stock Code: 870098



SAIKANGTM
LEADING BRAND OF CHINA
MEDICAL FURNITURE



LEADING BRAND OF CHINA
MEDICAL FURNITURE

RACE AGAINST TIME
HEALTHIER YOUR LIFE



CONTENTS

01	Strength of Saikang	09	ABS Series
03	Technology of Saikang	23	Phenolic Series
05	Innovation of Saikang	29	Color Painted/Stainless Steel Series
07	Concepts & Advantages	49	Accessories

STRENGTH OF SAIKANG

20

20 years focusing on design, production, service and technical consultation of medical care equipment

150

Located in Zhangjiagang City, 150km nearby Shanghai Economic Center

10000000

Saikang delivered 1,000,000 sets of products globally, Products all over Asia, Middle East, North America and Africa

400

There are more than 400 employees worldwide, more than 50 qualified engineers, and more than 100 trained technicians

189

Equiped with 189 sets of various automatic production facilities, Saikang follow strict standard control system on every production process, enable the traceable of each part

50000

The finished phase I production area covers more than 30,000 square meters, and phase II will be more than 20,000 square meters of advanced production area

16

Saikang offer comprehensive medical solutions with 16 product categories

562

Safe, efficient and easy-to-use medical solution including emergency, nursing and ICU wards of more than 500 series of products

29

With CE certification, ISO9001/13485 quality system certification, U.S. FDA and other 29 certifications

126

Ownership of 126 patents, continuously leading new standards in defining product quality and service efficiency

108

Since the establishment of the company, products and solutions have been used in more than 108 countries

1000

Saikang products are committed to creating a modern nursing environment, accumulatively serving more than 1,000 medical and health care organizations



LEADING BRAND OF MEDICAL FURNITURE

Patients in medical rehabilitation often have to face huge challenges. Many people who suffering from diseases or accidents are faced with retraining their minds and bodies and must adapt to the new environment. In this case, the positive atmosphere would offer decisive help. Hospital beds and furniture should provide the best support for patients and nurses. Considering this case, Saikang devoted to provide medical care solutions and perfect choices for medical institutions and rehabilitation centers, to help patients achieve rapid or long-term recovery.

Saikang was founded in 2002, now is a leading global provider of medical furniture and solutions. Head office locates in Suzhou, and branches in US, UK, Dubai, Morocco and Kenya, having 400 employees worldwide. Saikang has formed a large global network of research and development, marketing and service. Currently, Saikang products and medical solutions have been used in more than 108 countries and regions in the world.

TECHNOLOGY OF SAIKANG

SAIKANG insists on providing reliable quality, reshapes the "Made in China" products in worldwide market. SAIKANG has obtained ISO9001/13485 quality system certification, CE safety certification, EU safety certification, and US FDA market access license.

Alibaba.com Assessed Supplier

Supplier Assessment
Certificate Number: 779799_P-1
Presented to
Jiangsu Saikang Medical Equipment Co., Ltd.
江苏赛康医疗设备有限公司

For details, please refer to <http://saikang.en.alibaba.com>
Valid from: 06 May 2014 Valid until: 05 May 2015
Assessed by Mr. Marco Ma on 05 May 2014

CE Technical Documentation Review Report
Jiangsu Saikang Medical Equipment Co., Ltd.
Jiangsu Saikang Medical Equipment Co., Ltd.
Jiangsu Saikang Medical Equipment Co., Ltd.

Report Number: 1008462 001
Examination intent: Examination the completeness of the Technical Documentation according to the requirements of the Medical Devices Directive 90/269/EEC Annex VI

Product(s): SGP series, A series, SGT series
Type(s)/Model(s): SGP series, A series, SGT series
Classification: Class I, rule 1

Review result: During the examination of the provided Technical Documentation (the 3P/10/001, Ver.01, Date 2013-07-20), no non-compliance according to the requirements of the Medical Devices Directive 90/269/EEC Annex VI was detected.

Shanghai, 2015-08-13
TUV Rheinland (Shanghai) Co., Ltd.
Lead Auditor: Product Approval
Medical Device Services

Manufacturer: Jiangsu Saikang Medical Equipment Co., Ltd.
Report Number: 1008462 001
Examination intent: Examination the completeness of the Technical Documentation according to the requirements of the Medical Devices Directive 90/269/EEC Annex VI

Product(s): SGP series, A series, SGT series
Type(s)/Model(s): SGP series, A series, SGT series
Classification: Class I, rule 1

International Organization for Standardization



Patent



Eu CE



North American FDA



Министерство здравоохранения Российской Федерации



Saudi Food & Drug Authority



INNOVATION OF SAIKANG

With a total area of more than 30,000m², SAIKANG can fulfill various requirements of production. Through the intelligent management system, each process is visualized and standardized. The manufacturing base has its strict quality management system, so every production procedure is enabled to follow their standards and every produced piece is traceable.



Manufacturing Innovation

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality.




Maintenance Innovation

SAIKANG medical have more than 50 engineers, and qualified service term that composed of professional subcontractors, trained technicians, and after-sale staff. They are experienced and have excellent skills , not only to respond each service request on time, but also precaution every detail in equipment operation , always focus on products safely.



Service Innovation

SAIKANG medical have thorough service architecture: 7days x 24hours online customer call center, complete spare parts warehouse, advanced CRM system, precautionary service, etc. No matter where customer is, SAIKANG service will be always close at hand.



Design Innovation

SAIKANG medical designers have always maintain close communication with medical service providers, to help them providing better solutions to their patients. With long-standing experience in product structure and functional experience, SAIKANG continue to forge ahead with innovation and play a vital role in helping medical service providers to improve the quality of care.



Research Innovation

SAIKANG medical owns independent research and development center. Since its inception, every year SAIKANG invested 10% of its annual sales in R&D, launched 10 new products, and applied at least 10 patented technologies on an average basis, which lead SAIKANG to continuously define the quality of medical products and new standards of efficiency.

SAIKANG

LEADING BRAND OF CHINA
MEDICAL FURNITURE



CONCEPTS & ADVANTAGES



Standardization Certificates Approval

Obtain ISO9001, 13485, 14001, 18001 quality system certification, TUV approval for CE conformity, US FDA,GB national safety standard qualified enterprises.In accordance with the demand of facility management and environmental safety in JCI.



Ergonomic Design

As a member of the mobile terminal system, our medical trolley Have a high proportion of science, with a handiness easy push wheel design, make the handle more comfortable. The whole trolley adopts professional seamless design, easy to clean.



Nursing Intensity Optimal Design

In accordance with the demand of facility management and environmental safety in JCI. Each trolley with storage capacity can meet the demand of 8-10 beds per day.Trolleys with distinct signs to distinguish different using purposes, realize the classification management of non-medical waste and medical waste, which effectively improve the quality of medical care.





Chromatic Design

The design of medical carts integrated into the medical environment to regulate the mood of nurses and enhance the harmony between doctors and nurses. The chromatic design is according to the 2008 edition of "Guidelines for Hospital Construction in China" issued by the Hospital Management Institute of the National Health and Family Planning Commission,also according to "American Medical and Nursing Space" and "Health Nursing Art Program".

ABS SERIES

Leading Branding of China Medical Furniture
www.chinasaikang.com



ABS TROLLEY
SAIKANG MEDICAL

DAISY / SERIES

SKR054-ET
Emergency Trolley

Size: 850×520×960mm
750×475×930mm



SKR054-ET12
Emergency Trolley

Size: 850×520×960mm
750×475×930mm



SKR054-ET31
Emergency Trolley

Size: 850×520×960mm
750×475×930mm



DAISY / SERIES

SKR054-MT11
Medicine Trolley

Size: 850×520×960mm
750×475×930mm



SKR054-MT
Medicine Trolley

Size: 850×520×1010mm



SKR054-AT313
Anesthesia Trolley

Size: 750×475×930mm



SKR054-AT
Anesthesia Trolley

Size: 850×520×960mm
750×475×930mm



SKR054-AT12
Anesthesia Trolley

Size: 850×520×960mm
750×475×930mm



SKR054-AT21
Anesthesia Trolley

Size: 850×520×960mm
750×475×930mm



SKR054-AT212
Anesthesia Trolley

Size: 850×520×960mm
750×475×930mm



SKR054-CT91
Clinical Trolley

Size: 750×475×930mm



SKR054-CT91-1
Clinical Trolley

Size: 750×475×930mm



SKR054-IV
Infusion Trolley

Size: 750×475×930mm



SKR054-IV-1
Infusion Trolley

Size: 750×475×930mm



JASMINE / SERIES

SKR085-ET
Emergency Trolley

Size: 850×520×990mm



Main Body
ABS Plastic And Steel-plastic, Easy To Clean.




Table Top | Guard Rail | Push Handle


All-in-one table top, ABS table top with transparent dust-proof pad, convex design to avoid dropping down



ABS plastic drawers sliding smoothly and silently.
3*3 drawer dividers for arbitrary partitioning.
5 drawers with centralized lock and plug-in labels to avoid dust. Different colors options for drawer holder.



Two small drawers
panel height 80mm
inner size: 508×358×67mm



Two middle drawers
panel height 120mm
inner size: 508×358×108mm



One big drawer
panel height 238mm
inner size: 508×358×223mm

VIOLET / SERIES

SKR310-MT
Medicine Trolley

Size: 750×475×930mm
850×520×960mm



SKR310-ET
Emergency Trolley

Size: 750×475×930mm
850×520×960mm



SKR310-AT
Anesthesia Trolley

Size: 750×475×930mm
850×520×960mm



SKR310-CT
Clinical Trolley

Size: 750×475×930mm
850×520×960mm



SKR310-CT-1
Clinical Trolley

Size: 750×475×930mm
850×520×960mm



SKR310-IV
Infusion Trolley

Size: 750×475×930mm
850×520×960mm



SKR310-IV-1
Infusion Trolley

Size: 750×475×930mm
850×520×960mm



ROSEMARY / SERIES

SKR121-MT
Medicine Trolley

Size: 625×475×930mm
673×480×960mm



SKR121-ET
Emergency Trolley

Size: 625×475×930mm
673×480×960mm



SKR121-AT
Anesthesia Trolley

Size: 625×475×930mm
673×480×960mm



SKR121-CT
Clinical Trolley

Size: 625×475×930mm
673×480×960mm



SKR121-CT-1
Clinical Trolley

Size: 625×475×930mm
673×480×960mm



SKR121-IV
Infusion Trolley

Size: 625×475×930mm
673×480×960mm



SKR121-IV-1
Infusion Trolley

Size: 625×475×930mm
673×480×960mm



LILY / SERIES

SKR160-MT
Medicine Trolley

Size: 625×475×930mm



SKR160-ET
Emergency Trolley

Size: 625×475×930mm



SKR160-AT
Anesthesia Trolley

Size: 625×475×930mm



SKR160-CT
Clinical Trolley

Size: 625×475×930mm
673×478×930mm



SKR160-CT-1
Clinical Trolley

Size: 625×475×930mm
673×478×930mm



SKR160-IV
Infusion Trolley

Size: 625×475×930mm
673×480×960mm



SKR160-IV-1
Infusion Trolley

Size: 625×475×930mm
673×480×960mm



LOTUS / SERIES

SKR220-CP
Computer Cart

Size: 520×510×930mm



SKR220-CP31
Computer Cart

Size: 520×510×930mm



SKR220-CT-1
Clinical Trolley

Size: 1100×640×1100mm



SKR220-MT
Clinical Trolley

Size: 1100×640×1100mm



SKR220-CT81
Clinical Trolley

Size: 1100×640×1000mm



SKR220-CT
Clinical Trolley

Size: 1100×640×1000mm



VANILLA / SERIES

SKR023-7
Case History Trolley

Size: 450×410×960mm



SKR022-7
Case History Trolley

Size: 450×410×950mm



SKR022
Case History Trolley

Size: 450×410×950mm



SKR021-8
Case History Trolley

Size: 750×410×980mm



SKR021-7
Case History Trolley

Size: 750×410×980mm



SKR021-4
Case History Trolley

Size: 750×410×980mm



SKR-AB00
New All-in-One Computer Cart

Optional:
portable plug-in high-performance lithium battery
power supply system (with LED keyboard lights),
provide uninterruptible power for equipment or
computer, safe, reliable and convenient.

Casters:
Dual wheel casters for medical use, mute and
smooth, the front two casters equipped with double
pedal and total lock, 360 degree brake control, rear
followed by two free casters.

Base:
ABS plastic shell, stability, easy to clean and maintain,
high strength metal structure base frame, strong load
capacity and stability, the base integrated with foot
pedal to control height adjustment.



SKR-IB02
Simple Laptop Cart



SKR-J21
Endoscopic System
Equipment Cart

Size: 664×729×1666mm



SKR-A01
Simple All-in-One
Computer Cart



SKR-B00
Simple All-in-One
Computer Cart



SKR-J10
Patient Monitor Cart



SKR-J22
Patient Monitor Cart



MORNING GLORY / SERIES

SKR004
Stainless Steel-ABS Trolley
Size: 730×480×900mm



SKR004A
Stainless Steel-ABS Trolley
Size: 730×480×900mm



SKR004B
Stainless Steel-ABS Trolley
Size: 730×480×900mm



SKR005
Stainless Steel-ABS Trolley
Size: 730×480×900mm



SKR005A
Stainless Steel-ABS Trolley
Size: 730×480×900mm



SKR005B
Stainless Steel-ABS Trolley
Size: 730×480×900mm



IVY / SERIES

SKR001
ABS Utility Trolley
Size: 539×366×860mm
640×456×860mm
760×520×860mm



SKR001-1
ABS Utility Trolley
Size: 539×366×860mm
640×456×880mm
760×520×880mm



SKR002
ABS Utility Trolley
Size: 539×366×900mm
640×456×900mm
760×520×900mm



SKR002-1
ABS Utility Trolley
Size: 539×366×900mm
640×456×860mm
760×520×860mm



SKR001A
ABS Utility Trolley
Size: 539×366×860mm



SKR001A-1
ABS Utility Trolley
Size: 539×366×880mm



SKR002A
ABS Utility Trolley
Size: 539×366×900mm



SKR002A-1
ABS Utility Trolley
Size: 539×366×950mm



SKR001B
ABS Utility Trolley
Size: 640×456×860mm
760×520×860mm



SKR001B-1
ABS Utility Trolley
Size: 640×456×860mm



SKR002B
ABS Utility Trolley
Size: 640×456×900mm
760×520×900mm



SKR002B-1
ABS Utility Trolley
Size: 640×456×900mm
760×520×900mm






PHENOLIC SERIES

Leading Branding of China Medical Furniture
www.chinasaikang.com

Phenolic Resin Board



Structural Representation



Thickness Optional

Standard Thickness	Standard Thickness	Standard Thickness	Standard Thickness	Standard Thickness	Standard Thickness	Standard Thickness
mm	mm	Mm	mm	mm	mm	mm
4	6	8	10	12	13	16

Advantages

Space-saving

Colorful

Impact Resistance

Fire-proof

Water-proof

Chemical Resistance

GALAXY / SERIES

SKR701
Emergency Trolley

Size: 641×480×1121mm
641×480×902mm



SKR702
Emergency Trolley

Size: 588×480×971mm



SKR705
Infusion Trolley

Size: 588×480×891mm
588×480×1500mm



SKR706
Infusion Trolley

Size: 588×480×891mm
588×480×1500mm



SKR707
Infusion Trolley

Size: 588×480×891mm
588×480×1500mm



SKR703
Anesthesia Trolley

Size: 588×480×1063mm
588×480×1621mm



SKR704
Anesthesia Trolley

Size: 588×480×971mm
588×480×1529mm



SKR708
Medicine Trolley

Size: 588×480×1213mm



SKR709
Medicine Trolley

Size: 641×480×891mm



SKR710
Medicine Trolley

Size: 1206×480×1063mm



GALAXY / SERIES

SKR711

Clinical Trolley

Size: 1200×600×1213mm



SKR712

Clinical Trolley

Size: 588×480×891mm



SKR716

Clinical Trolley

Size: 588×480×891mm



SKR717

Nursing Trolley

Size: 588×480×891mm
350×408×700mm(储物袋)



SKR713

Fan-shaped Instrument Trolley

Size: 1438×580×891mm



SKR714

Clinical Trolley

Size: 1090×600×891mm



SKR715

Clinical Trolley

Size: 1090×600×891mm



SKR718

Case History Trolley

Size: 670×400×1038mm



SKR719

Case History Trolley

Size: 380×400×1025mm



SKR720

Case History Trolley

Size: 680×400×1025mm





COLOR PAINTED STAINLESS STEEL SERIES

Leading Branding of China Medical Furniture
www.chinasaikang.com



COLOR

Colorful Humanized Design

Adopting the influence of color on people's emotions, breaking the traditional color feeling of medical space and medical supplies.

Create a healing color combination that enhances vitality and aims for health.

To improve the patient's healing comfort, relieve the pressure of medical staff, and inject warmth into the medical space.



RHOMBUS / SERIES

SKR801
Movable Nurse Care Station (All in one)

Wireless Electronic Chemical Solutions

- High performance and high security battery system
- AC and DC power supply automatic conversion function
- Large capacity and low weight battery design reduces the work intensity of nursing staff
- LCD screen timely power display and alarm, one-touch switch
- One-piece machine-specific bracket, can be connected, install the market mainstream brand any model one machine, or a full set of hospital-type one machine to achieve wireless electronic work
- All-in-one bracket has lift mounting, rotation and tilt function
- Telescopic built-in cable management
- Table size: 620*470mm
- Work surface height: 930mm

Nursing Function Solution

- 8-10 beds daily infusion storage space
- Dozens of varieties, specifications of various types of sterile items, aseptic supplies, medicines, care appliances, classification
- Hand disinfectant, sharp sharp box, glove box, etc
- Various types of medical waste, domestic garbage sorting and recycling



Battery System

High performance and high security battery system
AC and DC power supply automatic conversion function
Large capacity and low weight battery design reduces the work intensity of nursing staff



One-piece Machine-specific Bracket

Can be connected, install the market mainstream brand any model one machine, or a full set of hospital-type one machine to achieve wireless electronic work
All-in-one bracket has lift mounting, rotation and tilt function



Drawer

Large, medium and small drawers are optional
Wide stainless steel handle for improved grip comfort
Ultra-quiet track is smooth



Telescopic Cable

Telescopic, built-in cable management



Lcd Screen Timely Power

Display and alarm, one-touch switch



Dust Bin

Flip cover automatically closed
Lightweight and beautiful, flexible assembly and disassembly



Sharp Container

Freely adjustable
Size can be choose



Metal Cover

Independent, large capacity



Color Painted Nursing Trolley



Concealed Side Pull Drawer

Silent and smooth rails, replaceable sub-workbench




Lateral Drawer Lock

Hidden on the side, avoiding the key blocking the first drawer from closing, avoiding the key colliding with it



Adjustable Medicine Division

The size of the drug can be adjusted according to the size of the drug placed



Grid Crash Pad

Soft material to avoid drug and drawer bump



Cpr Rescue Board

Easy to access, easy to use

Color Painted Nursing Trolley

(Optional) size for drawers(Small, middle, big)
(Optional) cabinet
(Optional) accessories

Corrosion resistant metal materials, color of coating could be customized.
resin rubber wheels, rotate 360°easily, noiseless and stable.
With dirt storage container that effectively separate medical waste from domestic waste
Multiple models, multiple parameters and multiple configurations

SKR804
Medicine Trolley
Size: 640×450×900/980mm



SKR853
Anesthesia Trolley
Size: 670×450×970/1500mm



SKR845
Emergency Trolley
Size: 670×450×940mm



SKR872
Clinical Trolley
Size: 670×450×940mm



SKR870
Medicine Trolley
Size: 640×450×890/970mm



SKR871
Medicine Trolley
Size: 670×450×940mm



SKR850
Emergency Trolley
Size: 670×450×950mm



SKR822
Clinical Trolley
Size: 600×400×850mm



SKR826
Clinical Trolley
Size: 710×620×990mm



SKR828
Clinical Trolley
Size: 640×450×880/980mm



SKR873
Clinical Trolley
Size: 690×450×840/970mm



SKR869
Clinical Trolley
Size: 640×450×880/960mm



Case History Trolley

(Optional) size and quantity for drawers
(Optional) number of case
Corrosion resistant metal materials, color and coating could be customized

With smooth and noiseless track
Size about the trolley could be customized
With lock bar to ensure patient information security

SKR875
Size: 705×410×940mm



SKR874
Size: 385×410×940mm



SKH013
Stainless Steel Case History Trolley
Size: 630×430×900mm



SKH014
Stainless Steel Case History Trolley
Size: 340×400×900mm



SKR023
Chart Holder
Size: 230×320mm
230×355mm
320×320mm



SKE020-1
Step Stool
Size: 405×230/170×340mm



SKE020-4D
Step Stool
Size: 405×230×220mm



SKE020-10
Step Stool
Size: 400×300×90mm



SKE020-8
Step Stool
Size: 496×342×500-120mm
496×342×600-159mm
496×342×750-200mm



SKE020-15
Step Stool
Size: 460×250×80mm
460×250×120mm
460×250×165mm



SKE017-3
Surgery Stool
Size: φ 350×460-600mm



SKE017-2
Surgery Stool
Size: φ 290×470-670mm



SKH006C
Instrument Trolley

Disassembly rail, easy to assemble
Movable connector is streamline designed with elegant appearance
Adjustable height of rails
Optional layers of rails



SKH003-2
Size: 460×460×860mm



Size: 670×450×900mm

SKH004
Instrument Trolley
Size: 670×450×900mm



SKH005
Fan-shaped Instrument Trolley
Size: 1550×460×900mm



SKH007-1
Instrument Trolley
Size: 1000×500×900mm



SKH030
Treatment Trolley
Size: 650×450×930mm



SKH038-1
Mayo Trolley
Size: 600×400×870-1100mm



SKH038-2
Mayo Trolley
Size: 500×350×900-1300mm



SKH001
Treatment Trolley
Size: 670×450×900mm



SKH002
Dressing Trolley
Size: 670×450×900mm



SKH006-5
Treatment Trolley
Size: 460×460×860



SKH002T
Dressing Trolley
Size: 670×450×900mm



SKH029
Dressing Trolley
Size: 670×450×900mm



SKH022
Treatment Trolley
Size: 600×400×850mm



SKH032
Treatment Trolley

Size: 670×450×1700mm



SKH028
Instrument Trolley

Size: 670×450×1000mm



SKH006-1
Instrument Trolley

Size: 670×450×900mm



SKH021
Emergency Trolley

Size: 600×400×850mm



SKH018
Anesthesia Trolley

Size: 600×400×850mm



SKH020
Emergency Trolley

Size: 600×400×850mm



SKH004-1
Instrument Trolley

Size: 670×450×900mm



SKH006-2
Instrument Trolley

Size: 670×450×900mm



SKH006-102
Instrument Trolley

Size: 670×450×1000mm



SKH019
Emergency Trolley

Size: 600×400×850mm



SKH012
Insulated Food Cart

Size: 670×450×1700mm



SK027-1
Dressing Trolley

The bar size and height was ergonomic design, operating more comfortable
Folding waste area,saving more space,easy to operate,
Quality canvas cloth,more durable.
Size is optional,bigger size of the waste area



Size: 655×555×940mm



Size: 1110×555×940mm

SKH027
Dressing Trolley

Size: 900×450×1050mm



SKH027-9
Dressing Trolley

Size: 900×450×1050mm



SKH027-2
Dressing Trolley

Size: 800×590×890mm

High quality tarpaulin,
Preventing sewage leaking,
ABS lid can separate the bad
smell and insects
Pedal system controls the
lid,can prevent touching the
trolley directly



SKH040-1
Dressing Trolley

Size: 800×600×840mm

Bigger size bag,easier to put
and take out
Quality canvas cloth,more
durable,lid is optional
Arc bump wheels can protect
the trolley



Two bags

Solid Linen Trolley



SKH040(1)

Size: 800×600×840mm



SKH040(2)

Size: φ 600×840mm

Stainless Steel Waste Basin

SKR041(5)

Size: φ 490mm×810~1250mm



SKH041(6)

Size: φ 490mm×810~1250mm

SKH034-5

Stainless Steel Movable Waste Bucket

Size: φ 400mm× φ 300mm×165mm×14L




Department Of Surgery (Medical Auxiliary Room)





- 1 Flat Opening Shelf
- 2 Manual Washing Sink
- 3 Toilet Waste Pool
- 4 Waste Collecting Table For Transferring
- 5 Single Washing Unit
- 6 Double Waste Trolley
- 7 Waste Trolley
- 8 Shoes Shelf
- 9 Mop Holder
- 10 Mop Pool

Central Supply Room



SKH036
Stainless Steel Washing Sink



可选配件



Cabient

High quality transparent glass

Optional layers

Glass size optional

High safety lock

Corrosion resistance, anti-aging handle

Solid and non-deformation

SKH051
Stainless Steel Cabient
Size: 900×400×1800mm



SKH052
Stainless Steel Cabient
Size: 900×400×1800mm



SKH054
Stainless Steel Cabient
Size: 900×500/300×1750mm



SKH055
Stainless Steel Cabient
Size: 1000×500/300×1750mm



SKH051A
Powder Coated Steel Cabient
Size: 900×400×1800mm



SKH052A
Powder Coated Steel Cabient
Size: 900×400×1800mm



SKH058-5A
Powder Coated Steel Cabient
Size: 900×400×1800mm



SKH057A
Powder Coated Steel Cabient
Size: 900×400×1800mm



Work Table

Use high quality corrosion resistant material

Customized size

Parameter customized

Disassembled type optional

SKH070
Work Table
Size: 1800×700×800mm
1200×600×800mm



SKH071
Work Table
Size: 1800×700×800mm
1200×600×800mm



SKH071-1
Work Table
Size: 1800×700×800mm
1200×600×800mm



SKH077
Work Table
Size: 1200×800/300×1600mm



SKH060
Rotary Medicine Shelf
Size: φ 920×1780mm



SKH062
Two-sided Adjustable Medicine Shelf
Size: 800×460×2000mm



SKH122-2
Delivery Trolley

Size: 650×420×1500mm

- High quality corrosion resistant material
- Large aperture basket
- Optional layers



SKH078-1
Storage Rack

Size: 1200×600×1800mm

- Disassembled type
- High quality corrosion material
- Optional layers
- Customized parameters



SKH078
Storage Rack

Size: 1200×600×1800mm



SKH088
Swab Stand

Size: 900×500×1500mm



Transport Trolley

- High quality corrosion resistant material
- Disassembled type optional
- Size can be customized
- Resin castors with brakes, moving easily

SKH039
Stainless Steel Trolley

Size: 750×500×950mm



SKH008
Stainless Steel Trolley

Size: 800×480×760mm



SKH039-2
Stainless Steel Trolley

Size: 750×500×950mm



SKH010-1
Stainless Steel Trolley

Size: 900×550×900mm



ACCESSORIES

Leading Branding of China Medical Furniture
www.chinasaikang.com



P001-Centralized Lock



P002-Password Lock



P003-Disposable Lock



P004-IV Pole With 2 Hooks



P005-IV Pole With 4 Hooks



P006-Double-side Caster



P007-Directional Castor



P008-Universal Covered Castor



P009-Single-side Castor



P010-Sharp Container



P011-Sharp Container



P012-Waste Bucket



P013-Two Waste Bucket



P014-Mini Dust Bin



P015-Flip Trash Can



P016-Double Trash Can With Rotating Lid



P017-Middle Size Waste Bucket



P018-Trash Can



P019-Trash Can



P020-Large Storage Box



P021-Plastic Medicine Tray




P022-ABS Drawer



P023-ABS Drawer



P024-48 Grids Medicine Tray



P025-Four-color Medicine Tray



P026-Medicine Tray



P027-Side Box



P028-Stainless Steel Utility Basket



P029-Stainless Steel Classification Basket




P030-Basket



P031-Narrow Medicine Tray



P032-File Basket



P033-File Basket



P034-Utility Container



P035-Power Outlet & Hooks



P036-Hand Soap Holder



P037-Label Box



P038-Anesthesia Box



P039-Oxygen Cylinder Holder



P040-Cpr Board




P041-Scissors Holder



P042-Gloves Holder



P043-Adjustable Medicine Tray



P044-Grid Crash Pad

NOTE

All the images, specifications and illustrations belong to Saikang's proprietary products. Any form of imitation or copy is a violation of the rights of our company. The actual color of products maybe have color difference as shown in catalogue. Products appearance maybe have difference if when technique upgraded.
Copyright © 2019 SAIKANG All Rights Reserved.