

24hrs: 0086-18021231901
www.chinasaikang.com

SAIKANG (USA Branch)
Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (East Africa Branch)
Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
P.O Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (China Headquarter)
No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (UK Branch)
Uk Healthward International Limited
Owner(s) Address DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (Mid East Branch)
Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City,Dubai UAE
TEL: +20 128 8962497

SAIKANG (Russia Branch)
454080, Russian Federation,
the city of Chelyabinsk,
Showroom, street: Truda 174,
TEL: +7 (351) 214-55-66,
214-88-11,
214-15-16



Factory Tour



Website



SAIKANG
Leading Brand of China Medical Furniture

www.chinasaikang.com

E20-1

SAIKANG®

LEADING BRAND OF CHINA
MEDICAL FURNITURE

JIANGSU SAIKANG MEDICAL EQUIPMENT CO., LTD.
www.chinasaikang.com

Stock Code: 870098



SAIKANG[®]
LEADING BRAND OF CHINA
MEDICAL FURNITURE



LEADING BRAND OF CHINA
MEDICAL FURNITURE

RACE AGAINST TIME
HEALTHIER YOUR LIFE



CONTENTS

01	Strength of Saikang	09	Hospital Bed	53	Mattress
03	Technology of Saikang	29	Homecare Bed	55	Bedside Table
05	Innovation of Saikang	31	Orthopedic Bed	56	Overbed Table
07	Quality of Saikang	33	Children Bed	57	Medical Chair
		37	Baby Crib	64	Stretcher
		38	Bathing Bed	67	Stretcher Trolley
		39	Examination Table	68	Patient Trolley
		42	Gynecological Table	73	Medical Trolley
		45	Obstetric Bed	79	Surgical Instrument
		47	Operating Table	80	Cabinet/Shelf

STRENGTH OF SAIKANG

20

20 years focusing on design, production, service and technical consultation of medical care equipment

150

Located in Zhangjiagang City, 150km nearby Shanghai Economic Center

1000000

Saikang delivered 1,000,000 sets of products globally, Products all over Asia, Middle East, North America and Africa

400

There are more than 400 employees worldwide, more than 50 qualified engineers, and more than 100 trained technicians

189

Equiped with 189 sets of various automatic production facilities, Saikang follow strict standard control system on every production process, enable the traceable of each part

50000

The finished phase I production area covers more than 30,000 square meters, and phase II will be more than 20,000 square meters of advanced production area

16

Saikang offer comprehensive medical solutions with 16 product categories

562

Safe, efficient and easy-to-use medical solution including emergency, nursing and ICU wards of more than 500 series of products

29

With CE certification, ISO9001/13485 quality system certification, U.S. FDA and other 29 certifications

126

Ownership of 126 patents, continuously leading new standards in defining product quality and service efficiency

108

Since the establishment of the company, products and solutions have been used in more than 108 countries

1000

Saikang products are committed to creating a modern nursing environment, accumulatively serving more than 1,000 medical and health care organizations




LEADING BRAND OF MEDICAL FURNITURE

Patients in medical rehabilitation often have to face huge challenges. Many people who suffering from diseases or accidents are faced with retraining their minds and bodies and must adapt to the new environment. In this case, the positive atmosphere would offer decisive help. Hospital beds and furniture should provide the best support for patients and nurses. Considering this case, Saikang devoted to provide medical care solutions and perfect choices for medical institutions and rehabilitation centers, to help patients achieve rapid or long-term recovery.

Saikang was founded in 2002, now is a leading global provider of medical furniture and solutions. Head office locates in Suzhou, and branches in US, UK, Dubai, Russia and Kenya, having 400 employees worldwide. Saikang has formed a large global network of research and development, marketing and service. Currently, Saikang products and medical solutions have been used in more than 108 countries and regions in the world.

TECHNOLOGY OF SAIKANG

SAIKANG insists on providing reliable quality, reshapes them or "Made in China" products in worldwide market. SAIKANG has obtained ISO9001:2015 quality system certification, CE approval, FCC conformity, EU safety certifications, and US FDA, CE, RoHS, PSE, Russian FDA market access license.





INNOVATION OF SAIKANG

Design Innovation



SAIKANG medical designers have always maintain close communication with medical service providers, to help them providing better solutions to their patients. With long-standing experience in product structure and functional experience, SAIKANG continue to forge ahead with innovation and play a vital role in helping medical service providers to improve the quality of care.

Research Innovation



SAIKANG medical owns independent research and development center. Since its inception, every year SAIKANG invested 10% of its annual sales in R&D, launched 10 new products, and applied at least 10 patented technologies on an average basis, which lead SAIKANG to continuously define the quality of medical products and new standards of efficiency.

Manufacturing Innovation



SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality.

Maintenance Innovation



SAIKANG medical have more than 50 engineers, and qualified service term that composed of professional subcontractors, trained technicians, and after-sale staff. They are experienced and have excellent skills , not only to respond each service request on time, but also precaution every detail in equipment operation , always focus on products safely.

Service Innovation



SAIKANG medical have thorough service architecture: 7days x 24hours online customer call center, complete spare parts warehouse, advanced CRM system, precautionary service, etc. No matter where customer is, SAIKANG service will be always close at hand.

QUALITY OF SAIKANG

Bed panels

A Harmonious Whole

Different panel design available, easy to remove and fix, without any tools. Fitted with buffers that protect both the walls and the bed itself when travelling.

Sleeping Surface


Unrivalled Sleeping Comfort

With 4 different sleeping surface to choose from. Profiled steel, thermoformed synthetic, shaped Plastic.

Protection Side Rails

Safety is Priceless

Atmosphere galleries, ABS collapsible side rails offer the patient protection on 2/3 of the bed. Thanks to the special design, the galleries constitute the ultimate protection for the patient.



Control Panel

Autonomous Control

Functional and easy to handle. The control offers all the functions which allow the patient and stuff to choose according the positions.

Height Adjustment System


Perfect Stability

The height-adjustable support arms guarantee excellent bed stability. Thanks to the special design [patented] the patient can more easily get out of bed and receive care.


World Class Branded Parts

Perfect Partners

We support and produce our products in accordance to EN international standard, our outsourced parts are from world class tracked companies.



Technique. Be Unique



Quality. Be Confident



Design. Be Harmonized



75°±10°
Back-rest

35°±10°
Knee-rest

450mm
Height
Adjustment
350-800mm

16°
Reverse
Trendelenburg

16°
Trendelenburg

Cardiac
Chair Position

Applicable
for C-arm
Examination

200mm
Bed Extension

weight
Weighing System

CPR
Instant Release

0-90°
Angle Indicator

100+60mm
Auto-regression

30°
Position Line

Built-in LED
Torch

Oasis Series

Y8y
Electric Bed
SKD-C

Technical Parameters

External Length	2130mm
External Width	1050mm
Bed Height Min.	350mm
Bed Height Max.	800mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Advantages

- ★ Double Auto-regression
- ★ Super-low Bed Position 350mm
- ★ PP Bed Platform, Removable & Washable
- ★ Inserted-type PP Head Foot Board, One Step Release With Safe Lock, Patent Design
- ★ PP Safety Guardrail With Two Steps Unlocking System, Patent Design
- ★ Back Rest 30degree Position, Reduce The Chance Of Lung Infection
- ★ Built-in LED Torch
- ★ Applicable For C-arm Examination

Technical Configuration

✓ Linear Motor	4pcs
✓ 5" Double Side Castors, Central Locking System	4pcs
✓ ABS Oasis Side Rail Extra Height	1set
✓ ABS Oasis Bed Ends with Safe Lock	1set
✓ PP Bed Platform, Removable & Washable	9pcs
✓ Plastic Socket Inside PP Guardrail for Phone or Glasses	2pcs
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ Ø100mm Wall Bumper Wheel	4pcs
✓ IV Pole	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ ABS Mattress Stopper	2pcs
✓ Hand Controller	1pc
✓ Manual CPR	

Optional

☐ Electric CPR
☐ Weighing System
☐ Bed Extension 20cm
☐ Back-up Battery
☐ Nurse Control Panel
☐ Emergency Stop Button
☐ Cardiac Chair Position
☐ Directional Wheels Function
☐ 5th and 6th Extra Side Rail For Full Protection
☐ Monkey Pole
☐ Balkan frame
☐ Linen Holder



360° View

Oasis Series

Y6y
Electric Bed
SKD-C

75°±10°
Back-rest

35°±10°
Knee-rest

440mm
Height
Adjustment
360-800mm

CPR
Instant Release


0-90°
Angle Indicator

200mm
Bed Extension

100+60mm
Auto-regression

30°
Position Line

Built-in LED
Torch



Technical Parameters

External Length	2130mm
External Width	1050mm
Bed Height Min.	360mm
Bed Height Max.	800mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Technical Configuration

☒ Linear Motor	3pcs
☒ 5" Double Side Castors, Central Locking System	4pcs
☒ ABS Oasis Side Rail Extra Height	1set
☒ ABS Oasis Bed Ends with Safe Lock	1set
☒ PP Bed Platform, Removable & Washable	9pcs
☒ Plastic Socket Inside PP Guardrail For Phone Or Glasses	2pcs
☒ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
☒ Ø100mm Wall Bumper Wheel	4pcs
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs
☒ Central Locking Pedal	2pcs
☒ ABS Mattress Stopper	2pcs
☒ Hand Controller	1pc

Optional

☐ Manual CPR	
☐ Electric CPR	
☐ Bed Extension 20cm	
☐ Back-up Battery	
☐ Nurse Control Panel	
☐ Emergency Stop Button	
☐ Directional Wheels Function	
☐ 5 th and 6 th Extra Side Rail For Full Protection	
☐ Monkey Pole	
☐ Balkan frame	

Oasis Series

Y6w
Electric Bed
SKD-C

75°±10°
Back-rest

35°±10°
Knee-rest

440mm
Height
Adjustment
365-805mm

200mm
Bed Extension

CPR
Instant Release

100+60mm
Auto-regression

Built-in LED
Torch



Technical Parameters

External Length	2130mm
External Width	1015mm
Bed Height Min.	365mm
Bed Height Max.	805mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Technical Configuration

☒ Linear Motor	3pcs
☒ 5" Covered Castors with Brake	4pcs
☒ Foldable Side Rail with ProErgo-Handle™	1set
☒ ABS Oasis Bed Ends with Safe Lock	1set
☒ PP Bed Platform, Removable & Washable	9pcs
☒ Ø100mm Wall Bumper Wheel	4pcs
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs
☒ ABS Mattress Stopper	2pcs
☒ Hand Controller	1pc

Optional

☐ Manual CPR	
☐ Electric CPR	
☐ Bed Extension 20cm	
☐ Back-up Battery	
☐ Nurse Control Panel	
☐ Emergency Stop Button	
☐ Directional Wheels Function	
☐ Monkey Pole	
☐ Balkan frame	



75°±10°
Back-rest

20°±10°
Knee-rest

450mm
Height Adjustment
500-950mm

18°
Reverse Trendelenburg

18°
Trendelenburg

weight
Weighing System

CPR
Instant Release

0-90°
Angle Indicator

100mm
Auto-regression

200mm
Bed Extension

Cardiac
Chair Position

X9x
Electric Bed
SKD-C

Technical Parameters

External Length	2040-2240mm
External Width	1020mm
Bed Height Min.	500mm
Bed Height Max.	950mm
Mattress Platform	1820-2020×820mm
Safe Working Load	250kg

Optional

- ☐ Weighing System With LCD Screen
- ☐ Data Storage System
- ☐ Monkey Pole

Technical Configuration

☒ Linak Motor	6pcs
☒ 6" Double Side Castors, Central Locking	4pcs
☒ ABS Side Rail with Integrated Control Panel	1set
☒ ABS Bed Ends One Step Removable	1set
☒ Ø75mm Wall Bumper Wheel	4pcs
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	4pcs
☒ 5" Directional Wheel Pedal	2pcs
☒ Central Locking Pedal	2pcs
☒ X-ray Translucent Mattress Platform	1set
☒ Back-up Battery	1set
☒ Night Light	1set
☒ Bed Push Handle for Emergency Transfer	1set
☒ Manual CPR	
☒ Electric CPR	



75°±10°
Back-rest

35°±10°
Knee-rest

290mm
Height
Adjustment
435-725mm

12°
Reverse
Trendelenburg

12°
Trendelenburg

CPR
Instant Release

0-90°
Angle Indicator

weight
Weighing System

Grace Care

D8d

Electric Bed
SKD-C

Technical Parameters

External Length	2280mm
External Width	1090mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

✓ Linear Motor	4pcs
✓ 5" Double Side Castors, Central Locking	4pcs
✓ ABS Grace Side Rail with Integrated Control Panel	1set
✓ ABS Grace Bed Ends with Integrated LCD Control Panel	1set
✓ Ø85mm Wall Bumper Wheel	4pcs
✓ Trapeze Socket	2pcs
✓ IV Pole	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ Manual CPR	

Optional

☐ Electrical CPR
☐ Weighing System with LCD Screen
☐ Bed Extension
☐ Back-up Battery
☐ Cardiac Position
☐ ABS Bed Surface
☐ X-ray Bed Surface
☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel

Grace Care

D6w

Electric Bed

SKD-C

75°±10°
Back-rest

35°±10°
Knee-rest

310mm
Height
Adjustment
435-745mm



Technical Parameters

External Length	2245mm
External Width	1000mm
Bed Height Min.	435mm
Bed Height Max.	745mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

☒ Linear Motor	3pcs
☒ Hand Controller	1pc
☒ 5" Hardshell Castors	4pcs
☒ Foldable Side Rail with ProErgo-Handle™	1set
☒ ABS Grace Bed Ends with Safe Lock	1set
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs

Optional

☐ Manual CPR
☐ Back-up Battery
☐ X-ray Bed Surface
☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel

Grace Care

D3d

Manual Bed

SK-CD1

75°±10°
Back-rest

35°±10°
Knee-rest

310mm
Height
Adjustment
435-745mm



Technical Parameters

External Length	2245mm
External Width	1090mm
Bed Height Min.	435mm
Bed Height Max.	745mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

☒ Manual Cranks	3pcs
☒ 5" Double Side Castors, Central Locking	4pcs
☒ ABS Grace Side Rail Extra Height	1set
☒ ABS Grace Bed Ends with Safe Lock	1set
☒ Central Locking Pedal	2pcs
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel

Grace Care

D3w

Manual Bed

SK-CD1

75°±10°
Back-rest

35°±10°
Knee-rest

310mm
Height
Adjustment
435-745mm



Technical Parameters

External Length	2245mm
External Width	1000mm
Bed Height Min.	435mm
Bed Height Max.	745mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

☒ Manual Cranks	3pcs
☒ 5" Hardshell Castors	4pcs
☒ Foldable Side Rail with ProErgo-Handle™	1set
☒ ABS Grace Bed Ends with Safe Lock	1set
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel

Grace Care

D2w

Manual Bed

SK-CD1

75°±10°
Back-rest

35°±10°
Knee-rest



Technical Parameters

External Length	2245mm
External Width	1000mm
Bed Height	510mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

☒ Manual Cranks	2pcs
☒ 5" Hardshell Castors	4pcs
☒ Foldable Side Rail with ProErgo-Handle™	1set
☒ ABS Grace Bed Ends with Safe Lock	1set
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs
☒ Shoes Holder	1set

Optional

☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel



75°±10°

Back-rest

35°±10°

Knee-rest

290mm

Height Adjustment

435-725mm

12°

Reverse Trendelenburg

12°

Trendelenburg

CPR

Instant Release

0-90°

Angle Indicator

B8e

Electric Bed

SKD-C

Technical Parameters

External Length	2170mm
External Width	1065mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

☒ Linear Motor	4pcs
☒ 5" Double Side Castors, Central Locking	4pcs
☒ ABS Side Rail Extra Height	1set
☒ ABS Bed Ends with Integrated Control Panel	1set
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs
☒ Central Locking Pedal	2pcs

Optional

☐ Manual CPR
☐ Electric CPR
☐ Weighing System
☐ Back-up Battery
☐ X-ray Bed Surface
☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel

B6e
Electric Bed
SKD-C


75°±10°
Back-rest

35°±10°
Knee-rest

310mm
Height
Adjustment
435-745mm

CPR
Instant Release

0-90°
Angle Indicator



Technical Parameters

External Length	2170mm
External Width	1060mm
Bed Height Min.	435mm
Bed Height Max.	745mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

☒ Linear Motor	3pcs
☒ Hand Controller	1pc
☒ 5" Hardshell Castors	4pcs
☒ ABS Side Rail Extra Height	1set
☒ ABS Bed Ends with Safe Lock	1set
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs

Optional

☐ Manual CPR
☐ Electric CPR
☐ Back-up Battery
☐ X-ray Bed Surface
☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel

B4w
Manual Bed
SK-CD1

75°±10°
Back-rest

35°±10°
Knee-rest

250mm
Height
Adjustment
460-710mm

12°
Reverse
Trendelenburg

12°
Trendelenburg



Technical Parameters

External Length	2180mm
External Width	960mm
Bed Height Min.	460mm
Bed Height Max.	710mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

☒ Manual Cranks	4pcs
☒ 5" Hardshell Castors	4pcs
☒ Foldable Side Rail with ProErgo-Handle™	1set
☒ ABS Bed Ends with Safe Lock	1set
☒ IV Pole	1set
☒ IV Pole Prevision	4pcs
☒ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel



A6k

Electric Bed

SK-CD1

Technical Parameters

External Length	2110mm
External Width	960mm
Bed Height Min.	400mm
Bed Height Max.	710mm
Mattress Platform	1925×820mm
Safe Working Load	135kg

Technical Configuration

✓ Linear Motor	3pcs
✓ 5" Covered Castors	4pcs
✓ Foldable Side Rail	1set
✓ ABS Bed Ends with Safe Lock	1set
✓ IV Pole	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	4pcs
✓ ABS Mattress Holder	2pcs

Optional

☐ Backup Battery
☐ Manual CPR
☐ Monkey Pole
☐ Directional Wheel

A3k
Manual Bed
SK-CD1

75°±10°
Back-rest

35°±10°
Knee-rest

310mm
Height
Adjustment
400-710mm



Technical Parameters

External Length	2110mm
External Width	960mm
Bed Height Min.	400mm
Bed Height Max.	710mm
Mattress Platform	1925×820mm
Safe Working Load	135kg

Technical Configuration

✓ Manual Crank	3pcs	✓ Drainage Hook	4pcs
✓ 5" Covered Castors	4pcs	✓ ABS Mattress Holder	2pcs
✓ Foldable Side Rail	1set		
✓ ABS Bed Ends with Safe Lock	1set		
✓ IV Pole	1set		
✓ IV Pole Prevision	4pcs		

A2k
Manual Bed
SK-CD1

75°±10°
Back-rest

35°±10°
Knee-rest



Technical Parameters

External Length	2110mm
External Width	960mm
Bed Height	490mm
Mattress Platform	1925×820mm
Safe Working Load	135kg

Technical Configuration

✓ Manual Crank	2pcs	✓ IV Pole Prevision	4pcs
✓ 5" Covered Castors	4pcs	✓ Drainage Hook	4pcs
✓ Foldable Side Rail	1set	✓ ABS Mattress Holder	2pcs
✓ ABS Bed Ends with Safe Lock	1set		
✓ IV Pole	1set		

A1k
Manual Bed
SK-CD1

75°±10°
Back-rest



Technical Parameters

External Length	2110mm
External Width	960mm
Bed Height	500mm
Mattress Platform	1890×830mm
Safe Working Load	135kg


Technical Configuration

✓ Manual Crank	1pcs	✓ IV Pole Prevision	4pcs
✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Foldable Side Rail	1set	✓ ABS Mattress Holder	2pcs
✓ ABS Bed Ends with Safe Lock	1set		
✓ IV Pole	1set		



360° View

H6k
Electric Homecare Bed
SKD-C



Technical Parameters

External Length	2015mm
External Width	1050mm
Bed Height Min.	360mm
Bed Height Max.	800mm
Mattress Platform	1925×880mm


Technical Configuration

✓ Linear Motor	3pcs	✓ IV Pole	1pc
✓ 5" Hardshell Castors	4pcs	✓ IV Pole Prevision	4pcs
✓ Aluminium Alloy Foldable Side Rail	1set	✓ Drainage Hook	2pcs
✓ Wood Bed Ends	1set	✓ ABS Mattress Holder	2pcs
✓ Wood Bed Side Panel	1set		

GB8e

Electric Orthopedic Traction Bed

SKD-C



Technical Parameters

External Length	2250mm
External Width	1060mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1935×900mm
Safe Working Load	250kg

Technical Configuration

Linear Motor	4pcs	Drainage Hook	2pcs
5" Double Side Castors, Central Locking	4pcs	Central Locking Pedal	2pcs
ABS Side Rail With Angle Indicator	1set	Aluminium Balkan Frame	1set
ABS Bed Ends with Safe Lock	1set		
IV Pole	1set		
IV Pole Prevision	4pcs		

GT3k

Manual Orthopedic Traction Bed

SKG-4C




Technical Parameters

External Length	2270mm
External Width	970mm
Bed Height	500mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

Manual Cranks	3pcs	IV Pole Prevision	4pcs
5" Double Side Castors, Central Locking	4pcs	Drainage Hook	2pcs
Central Locking Pedal	2pc	Aluminium Balkan Frame	1set
Aluminium alloy ABS Head and Foot Board With Safe Lock	1set		
Aluminium Alloy Foldable Side Rail	1set		
IV Pole	1set		

CX2x
Children Bed
SK-C2



Technical Parameters

External Length	1800mm
External Width	900mm
Bed Height	650mm
External Height	1260mm
Mattress Platform	1760×800mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole	1pc
✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Full Length Side Rail	1set	✓ Buffer Air Spring	4pcs
✓ Plastic Bed Ends with Cartoon Bear	1set		
✓ IV Pole Prevision	4pcs		

CR0q
Children Bed
SK-A2




Technical Parameters

External Length	1600mm
External Width	740mm
Bed Height	500mm
External Height	1040mm
Mattress Platform	1580×660mm

Technical Configuration

✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Full Length Side Rail	1set		
✓ Integrated Bed Ends	1set		
✓ IV Pole	1pc		
✓ IV Pole Prevision	4pcs		

CR2q
Children Bed
SK-A2



Technical Parameters

External Length	1600mm
External Width	700mm
Bed Height	500mm
External Height	1000mm
Mattress Platform	1580×660mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole Prevision	4pcs
✓ Full Length Side Rail	1set	✓ Drainage Hook	2pcs
✓ Intergrated Bed Ends	1set	✓ Shoes Holder	1pc
✓ Anti-slip Bed Feet Cover	4pcs		
✓ IV Pole	1pc		

CD3q
Children Bed
SK-CD1





Technical Parameters

External Length	2030mm
External Width	1045mm
Bed Height Min.	500mm
Bed Height Max.	770mm
Mattress Platform	1705×900mm

Technical Configuration

✓ Manual Cranks	3pcs	✓ IV Pole	1pc
✓ 5" Hardshell Castors	4pcs	✓ Drainage Hook	2pcs
✓ Aluminium Alloy Full Length Side Rail	1set		
✓ ABS Grace Bed Ends with Safe Lock	1set		
✓ IV Pole Prevision	4pcs		

CT8k
Electric Children Bed
SKD-C




Technical Parameters

External Length	1700mm
External Width	710mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1485×640mm

Technical Configuration

✓ Linear Motor	4pcs	✓ IV Pole	1pc
✓ 5" Double Side Castors, Central Locking	4pcs	✓ Drainage Hook	2pcs
✓ Aluminium Alloy Side Rail with Extra Protection	1set		
✓ Aluminium Alloy Bed Ends	1set		
✓ IV Pole Prevision	4pcs		

CT2k
Children Bed
SK-C2



Technical Parameters

External Length	1800mm
External Width	730mm
Bed Height	500mm
Mattress Platform	1580×660mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole	1set
✓ Aluminium Alloy Side Rail with Extra Protection	1set	✓ IV Pole Prevision	4pcs
✓ Aluminium Alloy Bed Ends	1set	✓ Drainage Hook	2pcs
✓ Anti-slip Bed Feet Cover	4pcs	✓ Shoes Holder	1pc

X01-1

Baby Crib
SK-AD3



Technical Parameters

External Length	875mm
External Width	525mm
Bed Height Min.	650mm
Bed Height Max.	900mm

Technical Configuration

☒ Hydraulic Pump	1 pc
☒ Bassinet	1 pc
☒ 3" Double Side Castors	4pcs

X01-5

Baby Crib
SK-B3



X01
Stainless Steel



Powder Coated Steel

Technical Parameters

External Length	810mm
External Width	520mm
Bed Height	990mm

Technical Configuration

☒ Bassinet	1 pc	☒ IV Pole	1 set
☒ 3" Covered Castors	4pcs	☒ Utility Basket	1 pc
☒ IV Pole Prevision	2pcs	☒ Adjustable Air Spring	1 pc

X20

Pediatric Exam Table
SKD-A

Technical Parameters

External Length	1800mm
External Width	650mm
Bed Height	830mm



Z0n

Electric Bath Bed
SKD-A



Z0y
Hydraulic



Technical Parameters

External Length	1950mm
External Width	780mm
Bed Height Min.	600mm
Bed Height Max.	830mm

Technical Configuration

☒ Linear Motor	1 pc
☒ 5" Covered Castors	4pcs
☒ Stainless Steel Side Rail	1 set
☒ Stainless Steel Bed Ends	1 set

X12

Electric Examination Table
SKD-A

230mm

Height Adjustment
500-730mm



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	500mm
Table Height Max.	730mm

Technical Configuration

☒ Linear Motor	1pc
☒ Anti-slip Bed Feet Cover	4pcs
☒ Removable Pillow	1pc

X14

Hydraulic Examination Table
SK-BD1

75°

Back-rest

400mm

Height Adjustment
500-900mm



Technical Parameters

External Length	1940mm
External Width	700mm
Table Height Min.	500mm
Table Height Max.	900mm

Technical Configuration

☒ Hydraulic Pump	1pc
☒ Adjustable Air Spring	2pcs
☒ 5" Covered Castors	4pcs
☒ Hydraulic pedal	2sets

X26

Electric Examination Table
SKD-A

75°

Back-rest

370mm

Height Adjustment
480-850mm



Technical Parameters

External Length	1810mm
External Width	680mm
Table Height Min.	480mm
Table Height Max.	850mm

Technical Configuration

☒ Pillow	1pc	☒ Hand Controller	1pc
☒ Motor	1pc	☒ Paper Roll Holder	1pc
☒ Mechanical Lever to Adjust Back Rest	1pc		
☒ 2" Casters	2pcs		

X27

Electric Examination Table
SKD-A

75°

Back-rest

30°

Waist Section up

440mm

Height Adjustment
530-970mm



Technical Parameters

External Length	1970mm
External Width	690mm
Table Height Min.	530mm
Table Height Max.	970mm

Technical Configuration

☒ Pillow	1pc	☒ Hand Controller	1pc
☒ Motor	2pcs	☒ Paper Roll Holder	1pc
☒ Gas spring	2pcs		
☒ 4" Casters	4pcs		

X07a
Examination Table
SK-A1



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

☒ Anti-slip Bed Feet Cover	4pcs

X09
Examination Table
SK-B1




Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

☒ Anti-slip Bed Feet Cover	4pcs
☒ Mechanical Lever to Adjust Back Rest	1set

A048
Gynecological Examination Table
SKF-B




Technical Parameters

External Length	1840mm
External Width	600mm
Table Height	800mm

Technical Configuration

☒ Drawer	4pcs	☒ Lithotomy Rod	1set
☒ Step Stool	1pc	☒ Stainless Steel Basin	1pc
☒ Retractable Foot Board	1pc	☒ Paper Roll Holder	1pc

A046
Gynecological Table
SKF-B




Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	800mm

Technical Configuration

☒ Crank Wheel	2pcs	☒ Stirrup	1set	☒ Anti-slip Bed Feet Cover	4pcs
☒ Lithotomy Rod	1set	☒ Hand Bar	1set	☒ Mechanical Lever	2pcs
☒ Shoulder Holder	1set	☒ Arm Holder	1set	☒ 3" Castors with Break	4pcs

A045-1
Gynecological Table
SKF-A




Technical Parameters

External Length	1820mm
External Width	500mm
Table Height	800mm

Technical Configuration

☒ Lithotomy Rod	1set	☒ Mechanical Lever	2sets
☒ Hand Bar	1set	☒ Stainless Steel Basin	1pc
☒ Anti-slip Bed Feet Cover	4pcs		

A99-8
Electric Gynecological Exam Couch
SKL-B



Technical Parameters

External Length	1450mm
External Width	600mm
Table Height	580mm

Technical Configuration

☒ Hand Bars	2pcs
☒ Leg Support with Cushion	2pcs
☒ Foot Pedal	1pc
☒ Stainless Steel Basin	1pc

A99-1
Electric Obstetric Bed
SKF-C



Technical Parameters

External Length	2000mm
External Width	600mm
Table Height Min.	650mm
Table Height Max.	900mm

Technical Configuration

☒ Bed frame	1pc	☒ Knee Curtch	2pcs
☒ Mattress	1pc	☒ Pedal Remote Controller	1pc
☒ Pull Handle	2pcs	☒ Stainless Steel Basin	1pc
☒ Armrest	2pcs		

A99-7

Electric Obstetric Bed

SKF-C

250mm

Height Adjustment

55°/18°

Back Plate up/down

35°/5°

Seat Plate up/down

90°

Leg Plate Swing Out



Technical Parameters

External Length	1980mm
External Width	720mm
Table Height Min.	630mm
Table Height Max.	880mm

Technical Configuration

☒ Bed Frame	1 pc	☒ Hand Bars	2 pcs	☒ Caster	1 pc
☒ Foot Pedal	1 pc	☒ Leg Support with Cushion	2 pcs		
☒ Hand Controller	1 pc	☒ Arm Support with Cushion	2 pcs		
☒ Anesthesia Screen	1 pc	☒ Stainless Steel Basin	1 pc		

A105

Manual Obstetric Table

SKF-C

200mm

Height Adjustment

50°/3°

Back Plate up/down

20°

Reverse Trendelenburg

10°

Trendelenburg

90°

Leg Plate Swing Out



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	700mm
Table Height Max.	900mm

Technical Configuration

☒ Bed Frame	1 pc	☒ Leg Support with Cushion	2 pcs
☒ Mattress	1 pc	☒ Arm Support with Cushion	2 pcs
☒ Anesthesia Screen Frame	1 pc	☒ Stainless Steel Basin	1 pc
☒ Hand Bars	2 pcs		

A98-3

Electric Obstetric Bed

SKL-A

70°

Back-rest

370mm

Height Adjustment

5°

Reverse Trendelenburg

12°

Trendelenburg



Technical Parameters

External Length	2050mm
External Width	980mm
Table Height Min.	580mm
Table Height Max.	950mm

Technical Configuration

☒ Leg Holder	1 pair	☒ IV pole	1 pc	☒ Castors	1 set
☒ Clamp	1 pair	☒ Controller	1 pc	☒ Stainless Steel	1 pc
☒ Grab Handle	1 pair	☒ Steel Guradrial	1 pair		
☒ Basin	1 pc	☒ Cushion	1 set		



A100-4A
Electric Operating Table
SKL-A

25°
Reverse
Trendelenburg

25°
Trendelenburg

15°
Lateral Tilt
left/right

40°/90°
Head Plate
up/down

80°/20°
Back Plate
up/down

15°/90°
Leg Plate
up/down

110mm
Kidney
Bridge Up

350mm
Horizontal
Sliding
C type Arm

Technical Parameters

External Length	2100mm
External Width	550mm
Table Height Min.	710mm
Table Height Max.	1010mm

Technical Configuration

Bed Frame	1 pc	Arm Holder with Cushion	2 pcs	Kidney Bridge Handle	1 pc
Mattress	1 set	Leg Support with Cushion	2 pcs		
Anesthesia Screen Frame	1 pc	Foot Support with Cushion	2 pcs		
Shoulder Holder with Cushion	2 pcs	Hand Controller	1 pc		



360° View

A100-4
Electric Operation Table
SKL-A



Technical Parameters

External Length	2070mm
External Width	550mm
Table Height Min.	700mm
Table Height Max.	1000mm

Technical Configuration

✓ Bed Frame	1pc	✓ Arm Holder with Cushion	2pcs
✓ Mattress	1set	✓ Leg Support with Cushion	2pcs
✓ Anesthesia Screen Frame	1pc	✓ Foot Support with Cushion	2pcs
✓ Shoulder Holder with Cushion	2pcs	✓ Kidney Bridge Handle	1pc



Low-position: 700mm



High-position: 1000mm



(Optional) Head Rest



Longitudinal Adjustment: 0-300mm



Longitudinal Adjustment: 0-300mm



(Optional) Lateral Arm Support



Reversed Trendelenburg: 0-25°



Trendelenburg: 0-25°



(Optional) Waste Basin For Gynecology



Kidney Bridge Elevation: 0-120mm



Anorectum Operation Position



E.N.T Operation Position



(Optional) Orthopedics Traction Frame



Lateral Tilt Right: 0-15°




Lateral Tilt Left: 0-15°

A106-2

Electro-Hydraulic Operation Table

SKL-A



Technical Parameters

External Length	2080mm
External Width	550mm
Table Height Min.	500mm
Table Height Max.	800mm


Technical Configuration

☒ Anaesthesia Screen	1 pc	☒ Foot Support	1 pc	☒ Fixing Clamp	7 pcs
☒ Hand Remote	1 pc	☒ Mattress	1 set	☒ Power Line	1 pc
☒ Kidney Bridge Handle	1 pc	☒ Shoulder Support	1 pair		
☒ Arm Support	1 pc	☒ Long Fixing Clamp	1 pair		

A3008A-1

Operation Table

SKT-C



Technical Parameters

External Length	2100mm
External Width	480mm
Table Height Min.	750mm
Table Height Max.	950mm

Technical Configuration

☒ Anaesthesia Screen	1 pc	☒ Foot Support	1 pc	☒ Fixing Clamp	8 pcs
☒ Hand Remote	1 pc	☒ Mattress	1 set	☒ Power Line	1 pc
☒ Kidney Bridge Handle	1 pc	☒ Shoulder Support	1 pair		
☒ Arm Support	1 pc	☒ Long Fixing Clamp	1 pair		

A03

Multi-Functional Orthopedic Frame

SKT-C


Technical Parameters

External Length	1700-1900mm
External Width	480mm
Table Height Min.	770mm
Table Height Max.	1020mm



Mattress (OEM Size)

SKP011



SKP104



SKP008



SKP006

Air Mattress



Mattress

Arrangement	2.5" bubble pad	Dimension	2000×900×70mm
Material	medical grade PVC	Weight	3kg

SKP013

Air Mattress



Mattress

Standard Size	80"×32"×5"	Bottom	tarpaulin /nylon PU
Top cover	Nylon PU /strength	Mattress weight	5.5-7kg
Air cells	Nylon / PVC/TPU	Weight Capacity	135-145kg

SKP012

Air Mattress



Mattress

Standard Size	80 "×32-36"	Bottom	tarpaulin /nylon PU
Top cover	Nylon PU /strength	Mattress weight	5.5-6.5kg
Air cells	Nylon / PVC/TPU	Weight Capacity	135-145kg

SKP015

ICU Air Mattress



Mattress

Standard Size	80"×32"×8"	Bottom	tarpaulin /nylon PU
Top cover	Nylon PU /strength	Mattress weight	10kg
Air cells	Nylon / PVC/TPU	Weight Capacity	180kg

SKS002
ABS Bedside Table
470×470×750mm(L×W×H)



SKS002-W
ABS Bedside Table
470×470×790mm(L×W×H)



SKS002-1
ABS Bedside Table
480×480×760mm(L×W×H)



SKS002-2
ABS Bedside Table
480×480×760mm(L×W×H)



SKS003
ABS Bedside Table
480×480×870mm(L×W×H)



SKS005
Steel Bedside Table
495×400×790mm(L×W×H)



SKS008
ABS+Steel Bedside Table
500×510×705mm(L×W×H)



SKS013-1
ABS+Steel Bedside Table
500×575×750mm(L×W×H)
480×460×760mm(L×W×H)



SKS009-2
Phenolic Resin Bedside Table
510×480×750mm(L×W×H)
Dining Table Max. Height: 950mm



SKS006-1
Stainless Steel Bedside Table
475×390×780mm(L×W×H)



SKS020
Solid Wood Bedside Table
500×430×650mm(L×W×H)



SKS019
Solid Wood Bedside Table
500×430×770mm(L×W×H)



SKS025
Phenolic Resin Bedside Table
475×485×810mm(L×W×H)

SKH046-11
Overbed Table

Technical Parameters
800×395×855-1155mm

Technical Configuration

- ✓ Air Spring for Height Adjustment
- ✓ Phenolic Resin Table Top



SKH046-2
Overbed Table



Technical Parameters
890×440×765-1050mm

Technical Configuration

- ✓ Air Spring for Height Adjustment
- ✓ ABS Table Top with Cup Socket



SKH046-14
Overbed Table



Technical Parameters
890×400×895-1195mm

Technical Configuration

- ✓ Air Spring for Height Adjustment
- ✓ ABS Table Top with Cup Socket and Push Handle



SKH042
Overbed Table



Technical Parameters
800×410×730-1130mm

Technical Configuration

- ✓ Air Spring for Height Adjustment
- ✓ MDF Table Top



SKH046-1
Telescopic Dining Board



Technical Parameters
870-1020×350×28mm

Technical Configuration

- ✓ Hook



SKE001

Attendant Chair

Technical Parameters

Seat Size	640×720×940mm
Bed Size	640×1900×420mm

Technical Configuration

☒ Arm	1 set
☒ 2" Castors	4pcs
☒ 3" Castors	2pcs
☒ Anti-slip Feet Cover	4pcs
☒ Pillow	1 pc



SKE001-3

Luxury Attendant Chair

Technical Parameters

Seat Size	710×760×950mm
Bed Size	710×1940×420mm

Technical Configuration

☒ Wood Arm Frame	1 set
☒ 2" Castors	2pcs
☒ Anti-slip Feet Cover	2pcs
☒ Pillow	1 pc



SKE001-4

Luxury Attendant Sofa

Technical Parameters

Seat Size	1350×900×650mm
Bed Size	1900×1350×500mm



SKE014-5

Luxury Attendant Sofa

Technical Parameters

Seat Size	900×900×650mm
Bed Size	1900×900×500mm



SKE-120B

Electric Two Function Dialysis Chair

0°-65°
Back Section
Adjustment

10°-90°
Leg Section
Adjustment

Technical Parameters

1830×580×560mm

Technical Configuration

- ☒ Hand Controller
- ☒ Power Line
- ☒ 3" Silent Caster
- ☐ Small Dining Table



SKE-120A

Electric Two Function Dialysis Chair

0°-70°
Back Section
Adjustment

0°-60°
Leg Section
Adjustment

Technical Parameters

1880×580×560mm

Technical Configuration

- ☒ Hand Controller
- ☒ Power Line
- ☒ 3" Silent Caster
- ☐ Small Dining Table



SKE-180

Electric Three Function Dialysis Chair

-8° - 70°
Back Section
Adjustment

0°-75°
Leg Section
Adjustment

550-750mm
Height
Adjustable

Technical Parameters

1800×580×(550-750)mm

Technical Configuration

- ☒ Hand Controller
- ☒ Power Line
- ☒ 3" Silent Caster
- ☒ With CPR
- ☐ Small Dining Table



SKE-170A

Electric Three Function Dialysis Chair

0°-70°
Back Section
Adjustment

20°-75°
Leg Section
Adjustment


540-640mm
Height
Adjustable

Technical Parameters

1950×580×(540-640)mm

Technical Configuration

- ☒ Hand Controller
- ☒ Power Line
- ☒ 3" Silent Caster
- ☐ Small Dining Table



SKE-100A
Manual Dialysis Chair

20° - 70°
Back Section
Adjustment

0° - 90°
Leg Section
Adjustment

Technical Parameters
1900×700×550mm

Technical Configuration

- ☒ Hand Controller
- ☒ 3" Silent Caster
- ☐ Small Dining Table



SKE087
Reclining Chair

Technical Parameters

Seat Size	820×850×1070mm
Bed Size	820×1600×450mm



SKE091
Blood Donation Chair

Functions

Back-rest	0°-70°
Leg Rest Down Fold	0°-90°

Technical Parameters

Seat Size	560×1030×1200mm
Bed Size	560×1650×470mm

Technical Configuration

☒ Arm	2pcs
☒ Anti-slip Feet Cover	4pcs



SKE090
Blood Donation Chair

Technical Parameters

Seat Size	650×800×450mm
-----------	---------------

Technical Configuration

☒ Arm	2pcs
☒ Anti-slip Feet Cover	4pcs



SKE004-1
Transfusion Chair

Technical Parameters

Seat Size	670×800×1050mm
-----------	----------------

Technical Configuration

☒ Anti-slip Feet Cover	4pcs
☒ Wooden Armrest	2pcs
☒ IV Pole Previsions	2pcs
☒ IV Pole	1 set



SKE005
Transfusion Chair

Technical Parameters

Seat Size	750×1000×1050mm
Bed Size	750×1550×450mm

Technical Configuration

☒ Anti-slip Feet Cover	4pcs
☒ Wooden Armrest	2pcs
☒ IV Pole Previsions	2pcs
☒ IV Pole	1 set
☒ Utility Basket	1 pc



Treat-waiting Chair



SKE008-1

Material: Powder Coated
1780×740×800mm(L×W×H)



SKE008

Material: Chromed Steel
1800×740×800mm(L×W×H)



SKE006-1

Material: PVC
1655×675×870mm(L×W×H)



SKE006-2

Material: PU
1765×675×905mm(L×W×H)





SKZ412
Office Table
Size: 1800×800×750mm



SKZ412-1
Side Cabinet
Size: 800×400×650mm



SKZ413
Meeting Table
Size: Ø800×750mm



SKH091
Cabinet
Size: 900×400×1800mm



SKE702
Office Chair
Mesh Cloth with Nylon Base
Size: 480×430×460-560mm



SKE055-8
Office Chair
Mesh Chair Back,high-quality Mesh with
Recycled Sponge, Plastic Handrail
Size: 480×470×900mm



SKE020-1
Step Stool
Size: 490×405×170/340mm



SKH049
Stainless Steel Screen
Size: 2000×1870mm



SK-CL001 **SK-CL501**
Curtain Track



SKH041-3
I.V. Stand
Size: 500×490×1310-2170mm



SKH041(1)
I.V. Stand
Size: Ø500×1350-2250mm



SKH034
Trash Can



12L

SKH034-5
Kick Bucket
Size: Ø400×Ø300×165×14L



14L



360° View

SKE180
Electric Wheelchair



SKE030
Wheelchair

Powder Coated Steel Structure
Seat Dimension: 395×455mm
Size: 1030×650×870mm



SKE031
Commode Chair

Height Adjustable, with Foot Pedal
Aluminium Folding Frame, Plastic
Size: 880×530×900mm



SKE053-2
Training Chair

Writing Board (370×245mm)
Size: 560×650×750mm



SKE052
Office Chair

Powder Coated Steel Structure
Back and Seat Section with the Elastic Fabrics
Size: 480×470×780mm



SKE711
Office Chair

Mesh Chair Back-high-quality Mesh with Recycled
Sponge, Plastic Handrail
Size: 500×490×500-600mm



SKE712
Office Chair

Mesh Chair Back-high-quality Mesh with Recycled
Sponge, plastic Handrail
Size: 500×400×1000mm



SKE802
Dental Chair

Chrome Plated Structure with Polyurethane Foam Seat
Size: Ø350×450-560mm



SKE015-2
Nurse Stool

Chrome Plated Structure with PU Foam Seat Top
Size: Ø340×490-680mm



SKE801
Nurse Stool

Stainless Steel Base with PU Leather,
High Density Sponge Seat
Size: Ø350×490-630mm



SKE015-5
Nurse Stool

Chrome Plated Structure with PU Leather,
High Density Sponge Seat
Size: Ø330×445-575mm



SKE013
Nurse Stool

Chrome Plated Structure with PU Foam Seat Top
Size: Ø330×445-575mm



SKE015-4
Nurse Stool

Chrome Plated Structure with PU Leather,
High Density Sponge Seat
Size: 440×400×450-590mmmm



SKE013-11
Lab Chair

Chrome Plated Structure with Polyurethane Foam Seat
Size: 330×230×450-590mmmm



SKE013-10
Lab Chair

Chrome Plated Structure with Polyurethane Foam Seat
Size: Ø350×435-575mm



SKE017-4
Stainless Steel Surgical Stool

201 Stainless Steel Structure with 430 Stainless
Steel Seat
Size: Ø295×500-700mm



SKB040(A001)
Foldable Stretcher
SKA-1C

Material: Aluminium Alloy Pipe & Canvas Fabric



Technical Parameters

Full Size	2100×500×100mm
Folded Size	580×220×210mm
Net Weight	6kg
Load Bearing	159kg



4-fold

SKB040(A002)
Foldable Stretcher
SKA-1B

Material: Aluminium Alloy Pipe & Canvas Fabric



Technical Parameters

Full Size	2200×530×140mm
Folded Size	1100×180×130mm
Net Weight	6kg
Load Bearing	159kg



2-fold

SKB040(A004)
Foldable Stretcher
SKA-1D

Material: Aluminium Alloy & Canvas Fabric



Technical Parameters

Full Size	1840×510×210mm
Folded Size	920×510×100mm
Net Weight	7.5kg
Load Bearing	159kg




(B) 4-fold



(A) 2-fold

SKB040(A006)
Scoop Stretcher

Material: Aluminium alloy



SKX-B1

SKX-B2

Technical Parameters

Full Size	2200×420×70mm
Folded Size	1660×420×70mm
Net Weight	8.5kg
Load Bearing	159kg



SPINE BOARD

SKB-1A



SKB2A10
1380×350×50mm



SKB2A14
1830×450×50mm



SKB2A11
1840×460×50cm



SKB2A12
1840×450×55mm



SKB2A13
1860×410×50mm



SKB2C01

Helicopter Stretcher
SKU-C1

Ideal for Lifting & Transporting by Helicopter



SKB2C01
SKU-C2

Technical Parameters

Full Size	2210×630×170mm
Load Bearing	272kg



SKB1C01-1

Foldable Transfer Stair Stretcher
SKW-B2




Technical Parameters

Full Size	850×520×1370mm
Folded Size	1100×540×180mm

SKB1C10

Electric Stair Stretcher
SKW-B2




Technical Parameters

Full Size	1230×500×1600mm
Folded Size	1020×550×210mm

SKB1C02-1

Manual Stair Stretcher for Downstairs
SKW-B2



Technical Parameters

Full Size	1250×550×1630mm
Folded Size	1250×550×330mm

Ambulance Stretcher Trolley

Material: Aluminium Alloy Foam Mattress: PVC Cover

SKB039(A)

SKC-2B1



Technical Parameters

Full Size	1900×550×400mm	Net Weight	25.5kg
Back-rest	70°	Load Bearing	159kg

SKB039(C)

SKC-3B4



Technical Parameters

Full Size	1900×550×920mm	Net Weight	36.8kg
Lowest Size	1900×550×260mm	Load Bearing	159kg
Back-rest	70°		

SKB039(E)

SKC-3F5



Technical Parameters

Full Size	1900×540×800mm	Net Weight	45kg
Lowest Size	1900×540×200mm	Load Bearing	159kg
Back-rest	85°		

SKB039(D)

SKC-3B5



Technical Parameters

Full Size	1950×550×870mm	Net Weight	40kg
Lowest Size	1950×550×230mm	Load Bearing	159kg
Back-rest	70°		

SKB037(A)

Stretcher Trolley

SKM-C

Size: 1960×600×810mm



SKB037(C)

Patient Trolley

SKM-G

Size: 1960×700×620-860mm



75°

Back-rest



35°

Knee-rest



240mm

Height

Adjustment

620-860mm



SKB038-4

Patient Trolley

SKM-G

Size: 1960×700×620-900mm



75°

Back-rest



280mm

Height

Adjustment

620-900mm



SKB037(B)

Patient Trolley

SKM-C

Size: 1960×600×770mm



75°

Back-rest



SKB038-2

Patient Trolley

SKM-E

Size: 1960×700×620-860mm



75°

Back-rest



240mm

Height

Adjustment

620-860mm



SKB041-3S

Hydraulic Patient Trolley

SKM-G

Size: 1950×650×620-950mm



55°

Back-rest



330mm

Height

Adjustment

620-950mm



SKB041-1

Patient Transportation Trolley

SKM-F

75°±10°
Back-rest

300mm
Height
Adjustment
540-840mm



Technical Parameters

Length	1930mm
Width	640mm
Min. Height	540mm
Max. Height	840mm

Technical Configuration

✓ Crank System	1 set	✓ Directional Wheel	1 pc	✓ Central Locking Break Pedal	2 pcs
✓ Back Rest Gas Spring	1 pc	✓ IV Pole Previsions	2 pcs	✓ Transfer Mattress	1 pc
✓ ABS Side Rail	1 set	✓ IV Pole	1 set	✓ Directional Pedal	2 pcs
✓ 6" Castors with Central Locking	4 pcs	✓ Oxygen Cylinder Holder	1 pc		

SKB041-2

Patient Transportation Trolley

SKM-F

75°±10°
Back-rest

300mm
Height
Adjustment
630-930mm

15°
Reverse
Trendelenburg

15°
Trendelenburg



Technical Parameters

Length	1930mm
Width	640mm
Min. Height	630mm
Max. Height	930mm

Technical Configuration

✓ Hydraulic Pump	2 pcs	✓ Directional wheel	1 pc	✓ Directional Pedal	2 pcs
✓ Back Rest Gas Spring	1 pc	✓ IV Pole Previsions	2 pcs	✓ Transfer Mattress	1 pc
✓ ABS Side Rail	1 set	✓ IV Pole	1 set		
✓ 8" Castors, Central Locking	4 pcs	✓ Central Locking Break Pedal	2 sets		

SKB041-3

Patient Transportation Trolley

SKM-F

75°±10°
Back-rest

45°
Knee-rest

300mm
Height
Adjustment
630-930mm

15°
Reverse
Trendelenburg

15°
Trendelenburg



Technical Parameters

Length	2100mm
Width	850mm
Min. Height	630mm
Max. Height	930mm

Technical Configuration

✓ USA Hydraulic Pump	2 pcs	✓ Ø74mm Wall Bumper Wheel	4 pcs
✓ Crank System	1 set	✓ IV Pole Previsions	4 pcs
✓ Back Rest Gas Spring	2 pcs	✓ Mattress	1 pc
✓ 8" Castors, Central Locking	4 pcs	✓ IV Pole	1 set
✓ Emergency Push Handle	1 set	✓ Central Locking Break Pedal	2 sets
✓ Directional Wheel	1 pc	✓ Directional Wheel Release Pedal	2 pcs
✓ Folding Side Rail	1 set	✓ X-ray Translucent Bed Platform	1 set

SKB041-6

Patient Transportation Trolley
SKM-F



Technical Parameters


Length	1930mm
Width	660mm
Min. Height	510mm
Max. Height	850mm

Technical Parameters

☒ Crank System	1set	☒ Directional Wheel	1pc	☒ Oxygen Cylinder Holder	1pc
☒ Back Rest Air Spring	1pc	☒ Directional Pedal	2pcs	☒ Central Locking Break Pedal	4pcs
☒ Plastic Side Rail	1set	☒ IV Pole Previsions	2pcs	☒ Transfer Mattress	1pc
☒ 6" Castors, Central Locking	4pcs	☒ IV Pole	1set	☒ Safe Lock for Side Rail	2pcs

SKB041

Operation Connecting Trolley
SKK-A



Technical Parameters

Length	3650mm
Width	640mm
Min. Height	640mm
Max. Height	970mm

Technical Configuration

☒ Crank System	2pcs	☒ 6" Castors with Central Locking	8pcs	☒ Directional Wheel Release Pedal	4pcs
☒ Back Rest Air Spring	1pc	☒ IV Pole Previsions	2pcs		
☒ ABS Side Rail	1set	☒ IV Pole	1set		
☒ Directional Wheel	2pcs	☒ Central Locking Break Pedal	4pcs		

SKR-AB00

New All-in-One Computer Cart



Worksurface	ABS Plastic Shell, High-strength Metal Frame
Function	Foot Pedal Control Height Adjustment
Castors	4" Double Side Castors
Optional	Portable Plug-in High-performance Lithium Battery

SKR-A01

Simple All-in-one Computer Cart



Worksurface	ABS Plastic Shell, Built-in Power Line
Frame	Applied For 22" All in One Computer
Castors	4" Double Side Castors

SKR-B00

Simple All-in-One Computer Cart



Worksurface	ABS Plastic Shell, Built-in Power Line
Frame	Applied For 22 Inch LCD Display
Castors	4" Double Side Castors
Optional	Built-in mini computer under worksurface

SKR-IB00

Simple Laptop Cart



Function	Handle Control Height Adjustment
Frame	Applied For ≤15" Laptop
Castors	4" Double Side Castors

SKR-J10

Patient Monitor Cart



Base	Integration Aluminum Alloy Column
Function	Fixed Height, Column with Attachments
Castors	3" Double Side Castors
Optional	Height can be Customized

SKR-J21

Endoscopic System Equipment Cart



Function	Endoscopic System Equipment Placing
Bracket	Height Adjustable, 360° Rotating
Layer	Four Layers, Two Middle Layer Can Be Adjusted
Castors	4" Double Side Casters

SKR054-ET

Emergency Trolley
Size: 750×475×920mm



Technical Configuration

✓ Defibrillator shelf	1 pc	✓ IV pole	1 pc	✓ Central Lock	1 set
✓ Sliding Side Shelf	1 pc	✓ CPR Board	1 pc	✓ File Box	1 pc
✓ Sharp Container	1 pc	✓ Oxygen Bottle Holder	1 pc		
✓ Dust Baskets	2pc	✓ Power Outlet	1 pc		

SKR051-ET

Emergency Trolley
Size: 750×475×930mm
850×520×960mm



SKR054-AT

Anesthesia Trolley
Size: 750×475×930mm
850×520×960mm



SKR054-MT

Medicine Trolley
Size: 750×480×920mm
850×520×950mm



SKR-CT750-1

Clinical Trolley
Size: 750×475×920mm



SKR-ET625

Emergency Trolley
Size: 625×475×930mm



SKR-AT625-1

Anesthesia Trolley
Size: 625×475×930mm



SKR-MT625

Medicine Trolley
Size: 625×475×930mm



SKR-CT625-1

Clinical Trolley
Size: 625×475×920mm



SKR004

ABS Trolley
Size: 730×480×900mm



SKR007

ABS Trolley
Size: 730×480×900mm



SKR021-8

Case History Trolley
Size: 750×410×950mm



SKR022

Case History Trolley
Size: 450×410×950mm



SKR001-2

ABS Trolley
Size: 539×366×845mm
640×456×845mm
760×520×845mm



SKR002-3

ABS Trolley
Size: 539×366×855mm



SK-MT650

Medicine Trolley
Size: 650×492×1040mm



SK-ET750

Emergency Trolley
Size: 750×530×1050mm



SK-ET75077A

Emergency Trolley
Size: 750×475×930mm



SK-AT75077A2

Anesthesia Trolley
Size: 750×475×930mm



SK-CT75077C2

Treatment Trolley
Size: 750×475×920mm



SKR-ET850

Emergency Trolley
Size: 850×520×990mm



SKR-CT600

Clinical Trolley
Size: 775×635×910mm



SKR-MT606

Clinical Trolley
Size: 775×630×1080mm



SKR-CT232

Clinical Trolley
Size: 780×475×1000mm



SKR-CT601

Clinical Trolley
Size: 775×635×910mm



SKR-CT670

Clinical Trolley
Size: 670×450×825-1000mm



SKR-CT670-1

Clinical Trolley
Size: 670×450×920-950mm



SKR-AT670

Anesthesia Trolley
Size: 670×450×880-960mm



SKR-ET670

Clinical Trolley
Size: 670×450×880-960mm



SKH040

Laundry Trolley
Size: 800×600×840mm



SKH040-1

Laundry Trolley
Size: 800×600×840mm



SKH027

Laundry Trolley
Size: 900×450×1050mm



SKR-WC561

Waste Collecting Trolley
Size: 800×590×890mm



SKH006C
Instrument Trolley
(Disassembly structure)
Size: 670×450×900mm



SKH006
Instrument Trolley
Size: 670×450×900mm



SKH001
Treatment Trolley
Size: 670×450×900mm



SKH002
Dressing Trolley
Size: 670×450×900mm



SKH006-1
Treatment Trolley
Size: 670×450×900mm



SKH004
Instrument Trolley
Size: 670×450×900mm



SKH006-5
Treatment Trolley
Size: 460×460×860mm



SKH005
Fan-shaped Instrument Trolley
Size: L: 1550×460×900mm
M: 1400×460×900mm
S: 1200×460×900mm



SKH029
Dressing Trolley
Size: 670×450×900mm



SKH018
Anesthesia Trolley
Size: 600×400×850mm



SKH032
Infusion Treatment Trolley
Size: 670×450×1700mm



SKH012
Insulated Food Cart
Heat store type
Size: 1320×580×1150mm



SKH021
Emergency Trolley
Size: 700×400×850mm



SKH019
Emergency Trolley
Size: 700×400×850mm



SKH028
Medicine Trolley
Size: 670×450×1000mm



SKH007-1
Stainless Steel Trolley
Size: L: 1000×500×900mm
M: 880×430×860mm
S: 760×360×820mm



SKH013
Stainless Steel Trolley
for Medical Record Holder (40 shelves)
Size: 630×430×900mm



SKH014
Stainless Steel Trolley
for Medical Record Holder (20 shelves)
Size: 340×400×900mm



SKH070
Work Table
Size: 1800×700×800mm



SKH077
Work Table
Size: 1200×800/300×1600mm



SKH078
Shelf With Four Layers
Size: 1200×600×1800mm



SKH105
Wet Linen Trolley
Size: 820×610×750mm



SKH048
Oxygen Bottle Trolley
Size: 1435×360×1040mm




SKH088
Swab Stand
Size: 900×500×1500mm




SKH036
Stainless Steel Washing Sink
Size: 1200(W)×600(D)×1100(H)mm



Size: 600(W)×600(D)×1100(H)mm



Size: 1800(W)×600(D)×1100(H)mm



Size: 2400(W)×600(D)×1100(H)mm



SKH041 (5)
Wash Basin
Size: Φ500×810-1300mm



SKH041 (6)

SKH038-2
Mayo Trolley
Size: 500×350×870-1300mm



SKN001

Stainless Steel Tank
Φ 170mm Φ 200mm Φ 230mm
Φ 260mm Φ 300mm



SKN002

Stainless Steel Cream Gauze Can
Φ 80×80mm Φ 90×90mm
Φ 100×100mm Φ 120×120mm
Φ 141×126mm Φ 160×154mm



SKN011

Stainless Steel Sterilization Tray
Without holes: 150×100×45mm
200×110×45mm
250×160×45mm
300×200×50mm
With holes: 150×100×45mm
200×110×45mm
250×160×45mm
300×200×50mm



SKN012

Stainless Steel Tray
300×200×48mm 350×250×48mm
400×300×48mm 450×300×48mm
350×500×48mm



SKN021

Stainless Steel Needle Sterilizer
Size: 77holes 105holes



SKN022

Stainless Steel Medicine-change Bowl
Φ 180mm Φ 200mm Φ 220mm
Φ 240mm Φ 260mm Φ 280mm



SKN031

Stainless Steel Kidney Dish
Shallow: L204×130mm
M190×116mm
S178×104mm
Deep: L243×125mm
M203×102mm
S165×80mm



SKN032

Stainless Steel Sterilized Bottle/
Medicine Cup
Sterilized bottles:
Φ 50×103mm Φ 50×93mm Φ 50×83mm Φ 50×53mm
Medicine cups: Φ 44×40mm Φ 50×50mm



SKN041

Stainless Steel Soak Forceps Bucket/
Measuring Cup/Enema Bucket
Measuring cup: 1L Φ 100×150mm
2L Φ 130×190mm
Enema bucket: 1L Φ 100×150mm
2L Φ 130×190mm



SKN042

Stainless Steel Forceps Bucket
Inverted cone: Φ 80×198mm
Φ 64×158mm
Φ 45×113mm
Straight: Φ 80×200mm
Φ 64×160mm
Φ 45×113mm



SKH090S

Cabinet
Size: 900×400×750mm



SKH090M

Cabinet
Size: 900×400×1050mm



SKH090L

Cabinet
Size: 900×400×1800mm



SKH092

Cabinet
Size: 900×400×1800mm



SKH093

Cabinet
Size: 900×400×1800mm



SKH098-2

Locker
Size: 900×500×1800mm



SKH098-3

Locker
Size: 900×500×1800mm



SKH098-4

Locker
Size: 900×500×1800mm




SKH098-6

Locker
Size: 900×500×1800mm




SKH098-9

Locker
Size: 900×500×1800mm



SKH059

Single-sided
Medicine Shelf
Size: 800×460×2000mm



SKH059-1
Size: 800×860×2000mm

SKH060

Rotary
Medicine Shelf
Size: 920×1780mm



SKH061

Adjustable Component
Medicine Shelf
Size: 1025×500×2000mm



SKH062

Two-sided Adjustable
Medicine Tray
Size: 800×860×2000mm



NOTE

All the images, specifications and illustrations belong to Saikang's proprietary products. Any form of imitation or copy is a violation of the rights of our company. The actual color of products maybe have color difference as shown in catalogue. Products appearance maybe have difference if when technique upgraded.
Copyright © 2020 SAIKANG All Rights Reserved.