



Stainless Steel Trolley

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

Addr.: No.35 Lehong Road, Modern Agriculture Demonstration Park,
Zhangjiagang City, Jiangsu Province, China.

Tel.: 0086-512-58525850

Fax.: 0086-512-58520713

P.C.: 215600

Web.: www.chinasaikang.com

E-mail.: inquiry@jssaikang.com





ABOUT US

SAIKANG MEDICAL is founded in 2002, combined with manufacture, sales, service, and support of some of the most innovative hospital furniture solutions. Managed by a professional team having a wide range of managerial experience, Saikang now is a leading provider of highly specialized medical furniture, such as electric and manual hospital beds, surgical tables, exam tables, infusion chairs, accompany chairs, hospital trolleys, stretcher and stretcher trolleys etc., registered in accordance with ISO9001:2000, ISO 13485:2003, CE, USFDA, Russia and Saudi registered and approved.



15 Years



24 Hour Service



400 Employees




Good Quality

Certificate



Our Workshop



Our Team



Stainless Steel Trolley

Treatment Trolley SKH006
Size: 670*450*900mm

Main Features:

- Stainless steel
- Two shelves
- Noiseless casters with brake
- Packing Size: 690*470*800(1 pcs /carton)
- Gross Weight: 12kg
- Net Weight: 11kg



Dressing Trolley SKH002
Size: 670*450*900mm



Treatment Trolley SKH001
Size: 670*450*900mm

Treatment Trolley SKH003-2



Dressing Trolley SKH002T
Size: 670*450*900mm

Instrument Trolley SKH004
Size: 670*450*900mm

Main Features:

- Stainless steel
- Three shelves
- With handle bar on short side
- Noiseless casters with brake
- Packing Method: Carton
- Gross Weight: 15kg
- Net Weight: 14kg



Instrument Trolley SKH004-1
Size: 670*450*900mm



Treatment Trolley SKH006-101
Size: 670*450*1000mm



Fan-shaped Instrument Trolley SKH005
Size: B: 1550*460*930mm M: 1400*460*900mm
S: 1200*460*900mm



Fan-shaped Instrument Trolley SKH005-1
Size: B: 1550*460*930mm M: 1400*460*900mm
S: 1200*460*900mm



Stainless Steel Trolley



Stainless Steel Trolley SKH007-1
Can buy it separately

- Main Features:**
- Stainless steel
 - With two shelves
 - 4 castors with brake



SKH007-1L



SKH007-1S



SKH007-1M



Treatment Trolley SKH006-1
Size: 670*450*900mm

- Main Features:**
- Stainless steel
 - Two shelves, with 2 drawers
 - Noiseless casters with brake
 - Packing Size: 690*470*800mm(1 pcs /carton)
 - Gross Weight: 17kg
 - Net Weight: 16kg



Treatment Trolley SKH006-5
Size: 46*46*86mm



Stainless Steel Trolley For Medical Record Holder (40 shelves) SKH013
Size: 630*430*900mm

- Main Features:**
- Stainless steel material
 - For holding 40 recorders
 - With top 2 drawers
 - With central lock
 - Noiseless casters with brake



Stainless Steel Trolley For Medical Record Holder (20 shelves) SKH014
Size: 340*400*900mm



- Main Features:**
- Stainless steel
 - Two shelves with two drawers
 - With one big drawer and cabinet with locks
 - Noiseless casters with brake

Anaesthetic Trolley SKH018
Size: 600*400*850mm



Emergency Trolley SKH021
Size: 600*400*850mm
With armrests: 700*400*850mm



Emergency Trolley SKH019
Size: 600*400*850mm
With armrests: 700*400*850mm

Stainless Steel Trolley



Dressing Trolley SKH029

Size: 670*450*900mm
With armrests: 860*450*900mm

Main Features:

- Stainless steel
- With two shelves
- Two drawers and tray
- Packing Size: 880*470*830(1 pcs /carton)
- Gross Weight: 14.5kg
- Net Weight: 14kg



Medicine Trolley SKH028

Size: 670*450*1000mm
With armrests: 790*450*1000mm



Infusion Treatment Trolley SKH032

Size: 670*450*1700mm



Wash Basin

SKH041(5) / SKH041(6)



Main Features:

- Stainless steel
- Tripod base
- A stainless steel sewage basin
- With 5 castors, with brake
- The size and height can be customized

Trash Can SKH034-5

Size: 14L



Main Features:

- All stainless steel structure
- With four castors, can movable

Oxygen Bottle Trolley SKH048

Size: 1435*360*1040mm



Stainless Steel Screen SKH049

Size: 2000*1870mm

Main Features:

- Four fold with curtains
- Stainless steel frame



SKH049-3
Three fold with curtains



SKH049-2
Two fold with curtains

Mayo Trolley SKH038-1

Size: 600*400*870-1100mm



Main Features:

- Stainless steel
- The top shelf with Tray its size: 600*400mm
- The height can be adjusted between 870mm and 1100mm

Mayo Trolley SKH038-2

Size: 500*350*900-1300mm



SKH037-1

Swab Stand SKH088

Size: 900*500*1500mm



Main Features:

- Size can be customized
- Stainless Steel Swab Stand with forty spikes i.e four transverse bars with 10 spikes each complete with stainless steel drip tray at base mounted on four antistatic castor wheels



Wet Linen Trolley SKH105

Size: frame: 765*555*420mm(highly wheeled)
Barrel diameter: 815*610*580mm

Laundry Trolley

Main Features:

- Anti-aging, anti-rust, erosion-resistant
- The bottom is equipped with silent wheel; easy to move and flexible
- Packing in a good quality thick paper carton(we can print your LOGO on the box)

Laundry Trolley SKH027
Size: 900*450*1050mm



Laundry Trolley SKH027-1
Size: 988*480*936mm



Main Features:

- Stainless steel
- Four casters, two layers with one drawer, handrail



Laundry Trolley SKH040
Size: 800*600*840mm

Laundry Trolley SKH040-1
Size: 800*600*840mm



Main Features:

- Stainless steel structure
- 4 antistatic castors with brake and bumpers
- With cavans bag, one white and one pink

Waste Collecting Trolley SKR-WC561
Size: 800*590*890mm



Main Features:

- Steel frame with plastic coated
- Big bag for waster, include a buffer system
- Noiseless operation
- With ABS plastic lid



Waste Collecting Trolley SKR-WC562
Size: 800*590*890mm

Work Table

Work Table SKH070
Size(L*W*H): 1200*600*800mm



Main Features:

- Stainless steel
- Stable Structure
- Easy for cleaning
- Multy use
- Size can be customized



SKH070C
Removable

Main Features:

- Stainless steel
- Stable Structure
- Easy for cleaning
- Multy use
- Size can be customized



Work Table SKH070-1
Size(L*W*H): 1800*700*800mm



SKH070C-1
Removable

Work Table SKH077
Size(L*W*H): 1200*800/300*1600mm



Main Features:

- Stainless steel
- Steel plate thickness: 0.8mm
- Stable table foot
- Thick material
- Size can be customized



Shelf With Four Layers SKH078
Size(L*W*H): 1200*600*1800mm

Step Stool

Step Stool SKE020-2

Main Features:

- Stainless steel frame
- Rubber step surface with stainless steel cover in case of skip
- Each layer height: 200 mm



Step Stool SKE020-2D

Main Features:

- Stainless steel frame
- Rubber step surface with stainless steel cover in case of skip



Step Stool SKE020-3

Main Features:

- Rubber step surface with stainless steel frame
- Plastic steel frame
- Each layer height: 200 mm



Step Stool SKE020-3D

Main Features:

- Rubber step surface with stainless steel frame
- Plastic steel frame



Step Stool SKE020-1

Main Features:

- Stainless steel frame
- ABS anti-skid stool surface(blue, gray)
- Size: 405*230/170*340mm



Step Stool SKE020-1D

Main Features:

- Stainless steel frame
- ABS anti-skid stool surface(blue, gray)
- Single steps
- Size: 405*230*220mm



Step Stool SKE020-4

Main Features:

- Plastic steel frame
- ABS anti-skid stool surface(blue, gray)
- Size: 405*230/170*340mm



Step Stool SKE020-4D

Main Features:

- Plastic steel frame
- ABS anti-skid stool surface(blue, gray)
- Size: 405*230*220mm

