



SK-L005

Examining Lamp



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE

Stock Code: 870098

TECHNICAL PARAMETERS

- / Arm: 105 +/- cm articulated, spring loaded arm, arm with on/off switch and incorporated electronically transformer
- / Mobile stand with swivel casters
- / Lamp emits natural white light: colour temperature approx. 4000 K
- / Halogen bulb approx. 12V/20W
- / Light intensity: approx 20.000 Lux at 40 cm
- / Reflector adjustable for positioning
- / Free cord: length approx 3m
- / Supplied with:
 - 1x spare bulb
 - 1x spare fuse
 - 1x Instructions for use



QUALITY CONTROL

ADJUSTABLE DOME

The lamp head angle is easy to adjust and the angle of adjustment is wide.



ADJUSTABLE HEIGHT

Height is adjustable, can match the height of the patient easily



POWER

The power cord is long enough



REFLECTOR

Reflector adjustable positioning



BASE

With swivel casters, easy to move





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

UK Healthward International Limited
Owner(s) Address:DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.