



# CR2q

Children Bed


SK-A2



**SAIKANG** / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE

Stock Code: 870098



Back-rest

Knee-rest

## TECHNICAL PARAMETERS

External L*W*H	1600*700*1000mm
Bed Height	500mm
Mattress Platform	1580*660mm
Back-rest Adjustment	0-75°(±10°)
Knee-rest Adjustment	0-35°(±10°)

## TECHNICAL CONFIGURATION

Full Length Side Rail	1 set
Intergrated Bed Ends	1 set
IV Pole	1 pc
IV Pole Prevision	4 pcs
Drainage Hook	2 pcs
Crank Handle	2 pcs

## STANDARD ACCESSORY

One telescopic IV pole

## OPTIONAL ACCESSORY

Basic foam mattress - 12cm thickness  
Anti-decubitus foam mattress - 14cm thickness  
Single arm lifting pole  
Utility Shelf

## QUALITY CONTROL



## RICH PARTS

LEADING BRAND OF  
CHINA MEDICAL FURNITURE




Punching platform, smooth surface no welding scar, bearing 250KG load. Option resin X-ray translucent platform



Powder coated steel frame, MDF material fire-proof plate. size of the head&foot board can be customized. With carton photograph



Foldable guardrail with protection handle, modern and beautiful, easy to install, easy to clean, anti-pinch, protect patient's safety



Plastic drainage hook, movable



Full SS IV pole, bearing 15kg



Shoes holder: size adjustable  
Oxygen cylinder rack - adjustable for anysize tank  
Lifting pole - Bearing 765N force  
Blocking







## INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

### JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901   
[www.chinasaikang.com](http://www.chinasaikang.com)

#### SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,  
Jiangsu Province,China.  
TEL: +86-512-58525850  
FAX: +86-512-58520713  
E-mail: [inquiry@jssaikang.com](mailto:inquiry@jssaikang.com)

#### SAIKANG (USA Branch)

Saikang Medical Technology, Inc  
17800 Castleton St Ste 180,  
City of Industry, Ca 91748  
Tel: +1 949648-2776

#### SAIKANG (UK Branch)


Uk Healthward International Limited  
Owner(s) Address DEPT 114,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +44 2032399738

#### SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.  
Viraj Complex Unit 3 Mombasa Road,  
PO Box 14262-00800, Nairobi.  
TEL: +254 733-704 690

#### SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC  
X22-S03-England Cluster,  
International City, Dubai UAE  
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.  
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.