




X12


Electric Examination Table

SKD-A



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098



Height Adjustment

TECHNICAL PARAMETERS

External L*W*H 1900*600* (500-800) mm

TECHNICAL CONFIGURATION

Linear Motor	1 pc
Anti-slip Bed Feet Cover	4 pcs
Removable Pillow	1 pc

ELECTRONICS

Voltage	100-240V
Frequency	50/60Hz
Plug	as required (hospital grade)



QUALITY CONTROL

LEATHER

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance

100,000

SPONGE

Flame retardant, 4 cm thickness, high-density foam, density 30 kg/m³

4. | 30

PAINTING

11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact

0.12. | 60° | 50.

METAL MATERIAL

Baosteel from Fortune Global 500
Panasonic Robotic ensure 360° full smooth welding

500 | 360°

RICH PARTS



1.2-1.5 mm Thickness powder coated steel frame, 250 kg bear loading



4 pcs Anti-slip bed feet covers



Hand controller, keep qualified after 2000 times use





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.