



SKB041


Operation
Connecting
Trolley



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE

Stock Code: 870098



Back-rest

Height Adjustment

MAIN TECHNICAL PARAMETERS

Size	3650*640*(640-970)mm
Back-rest	0-75°(±10°)

TECHNICAL CONFIGURATION

Crank System	2pc
Back Rest Air Spring	1pc
ABS Side Rail	1set
Directional Wheel	2pc
6" Castors with Central Locking	4pcs
IV Previsions	2pc
IV Pole	1set
Central Locking Break Pedal	4pcs
Directional Wheel Release Pedal	4pcs


PROTECT

cover all metal screws with rubber cap, protect screws from possible corrosion



CENTER WHEEL

Can be engaged to facilitate direction control, enabling users to navigate the stretcher easily in hectic and dynamic environments.



CONNECT WITH OPERATING ROOM

Transfer patient to operation room, prevent cross infection



QUALITY CONTROL

METAL MATERIAL

Baosteel from Fortune Global 500
Panasonic Robotic ensure 360°
smooth wedding



PLASTIC

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa



STEELS | PAINTING

1 | Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact



METAL DUST COVER

We use metal dust protection for crank system, more durable than plastic cover



BEARING

We use robust tube in weight bearing structure



AUTOLOCK

The patient stretcher will be autolocked when it is fixed to one side of the trolley



RICH PARTS

LEADING BRAND OF
CHINA MEDICAL FURNITURE



Adjustable vertical oxygen bottle, load capacity: 15KG



Stainless steel handle, more durable



Four double side casters of 6 inch diameter; simultaneously lockable by the step pedal at both sides.


TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts



Full SS IV pole, bearing 15kg



Water proof mattress with 2 belts





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.