



SKE008-1

TREAT-WAITING
CHAIR



SAIKANG / Medical


Healthier
your life

SKE008-1

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 1780×740×800 mm



Material

Powder coated steel.
Shining baked paint plating structure,
made of high quality cold rolled steel
plates.



Armrest

Made of high quality tensile steel
plate, with electroplated chromium
after derusting processing.



Beam

Made of high quality square tube, with
painting.



Backrest

Backrest with tilt, feel comfortable.



Color

Silver, color can be customized.



Optional

Optional PU seat, make people feel
more comfortable.





CONSTANT
INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE
BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue