




SKB2B01

SCOOP STRETCHER

Registration code: SKX-B1 / SKX-B2



SKX-B1



SKX-B2

SKB2B01


FEATURES AND HIGHLIGHTS

Technical parameters

▪ Product size (LxWxH)	2140x420x70mm
▪ Folded size (LxWxH)	1670x420x70mm
▪ Loading bearing	159kg
▪ Net weight	8.5kg

Design

It's a seperatable-type emergency stretcher for ambulance transfer fractured & serious patients, could strengthen the patients at origin place to reduce the further hurt to the patients.



Design

Equipped with a clutch hinge in the middle of both end, which can separate the stretcher into two parts




Material

Made of high-strength aluminum alloy materials



Length

The length of the scoop stretcher can be adjusted according to patients' body





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue