

X01-1

BABY CRIB

Registration code: SK-AD3



X01-1

FEATURES AND HIGHLIGHTS

Technical parameters

▪ External Length	875 mm
▪ External Width	525 mm
▪ Bed Height Min.	650 mm
▪ Bed Height Max.	900 mm

Technical configuration

<input checked="" type="checkbox"/> Hydraulic Pump	1pc
<input checked="" type="checkbox"/> Bassinet	1pc
<input checked="" type="checkbox"/> 3" Double Side Castors	4pcs

Design

Stylish modern design, ergonomic to provide a comfortable experience for babies.



Baby Basin

Transparent baby basin, environmentally plastic non-toxic, odorless, can be sterilized at temperatures above 100 °C.



Wheel

Rubber free wheel, wear-resistant, antiskid, shock absorption, can rotate 360° freely.



Handle

Easy to operate handle, control gas spring by handle to achieve various functions, excellent design and easy operation.



Color

Available in three colors: Blue, Green, Pink.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue