



A99-12

ELECTRIC OBSTETRIC BED

Registration code: SKF-C



SAIKANG / Medical

Healthier
your life

A99-12

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External (LxW):	1800x600mm
▪ Table height:	700-950mm
▪ Back plate up folding:	0-60°
▪ Reverse Trendelenburg:	0-25°
▪ Trendelenburg:	0-16°
▪ Load bearing:	≤150kg



Technical configuration

<input checked="" type="checkbox"/> Arm support	2 pcs
<input checked="" type="checkbox"/> Handle	2 pcs
<input checked="" type="checkbox"/> Leg support	2 pcs
<input checked="" type="checkbox"/> Stirrup	2 pcs
<input checked="" type="checkbox"/> Slider	4 pcs
<input checked="" type="checkbox"/> Tightener	2 pcs
<input checked="" type="checkbox"/> Hand controller	1 pc
<input checked="" type="checkbox"/> Foot switch	1 pc
<input checked="" type="checkbox"/> Stainless steel basin	1 pc



Material

The chassis cover is made of high quality 304 stainless steel, which is welded by the robot seamlessly. The shape is beautiful and stable.

**Mattress**

Polyurethane foam mattress, waterproof and seamless, easy cleaning and disinfection.

**Hand controller**

All the functions are controlled by hand controller.

**Braking system**

The pedal controls the brake function of the casters.

**Control System**

Electric control the height lifting, trendelenburg and reverse trendelenburg, back plate folding up, easy and convenient operation.

**Battery**

The delivery bed is equipped with a standard battery, which could support the electric operation for 50-80 times, and could keep maintaining the operation of various functions required by the operator under special circumstances such as power failure.

**Pull handles**

Both sides of the delivery bed are equipped with collapsible handle, which is convenient for the puerpera to get on and off the bed, with anti-slip lines, to prevent slipping when the force is applied.

**Sewage basin**

The delivery bed is equipped with a stainless steel basin, which can be stored under the bed through a high-quality three-stage slide rails.

**Auxiliary Table**

The auxiliary table can be retractable and stored under the bed body, easy to use.

**Leg supports**

The leg support is integrated with stirrup function, and the delivery bed is equipped with a pair of stirrups.

**Arm Support**

The arm support could be adjusted into multiple angles, it can be used to comfortably fix the puerpera's arm at various angles.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue