



# R000w

## ONE FUNCTION MANUAL BED

Registration code: SK-CD1




**SAIKANG** Medical

Healthier  
your life

# R000w

## FEATURES AND HIGHLIGHTS



### Technical parameters

▪ External size (LxWxH)	2000x1000x490 mm
▪ Mattress platform	1950x950 mm
▪ Safe working load	200 kg
▪ Back-rest adjustment	0-75° (±10°)



### Technical configuration

<input checked="" type="checkbox"/> Powder coated Bed Ends	1set
<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> IV Pole Prevision	1p

### Optional

<input type="checkbox"/> 5" Noiseless Castors	
<input type="checkbox"/> Integrated detachable guardrail	



**Steels | Painting**


Baosteel from Fortune Global500.  
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12 mm, brightness 60°, paint can resist 50kg impact.

**Welding**

Panasonic robotic ensure, 360° full smooth welding

**Bed Surface**

Reticulated bed platform.

**Head & Foot Board**

Powder coated head&foot board, simple installation,Removable.

**Handle**


Equipped with a handle to manually control the angle of back lift.

**IV Pole**

Full stainless steel IV pole, bearing 15 kg.

**Optional**

Integrated detachable guardrail.

**Optional**

5" noiseless castors.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue