



# SKB039(D)

## AMBULANCE STRETCHER TROLLEY

Registration code: SKC-3B5



# SKB039(D)

## FEATURES AND HIGHLIGHTS



### Technical parameters

▪ Full size (LxWxH)	1900x550x920mm
▪ Lowest size (LxWxH)	1900x550x250mm
▪ Back-rest lifting	0-75°
▪ Net weight	35kg
▪ Load bearing	≤159kg

#### Safety Lock

The safety lock is locked when transporting a patient, and unlocked when getting on an ambulance.



#### Caster

5" silicone caster. Advanced protection universal brake wheel, easy to move, simple to operate, reliable performance.



#### IV Pole

Telescopic stainless steel infusion rod, strong bearing capacity, more convenient to use.



#### Function

The mechanical folding legs are controlled by the handles at both ends of the stretcher, which can be easily lowered when entering the ambulance.



#### Back Rest Lifting

Back is adjustable which makes the wounded feel comfortable. Back rest lifting 0-75°.



#### Fixing Devices

Can be locked by the fixing devices when on the ambulance car.



#### Structure

Adopt high-quality aluminum alloy material, fitting aluminum alloy die-casting, light weight, anti-corrosion and anti-rust.



#### Cushion

Adopt sponge cushion to make the injured and sick lying comfortably.



#### Leather / Safety Belts

Artificial leather, bright appearance, soft texture, abrasion resistance, folding resistance, acid and alkali resistance. With 2 safety belts, make patient safer.



#### Guardrail

The aluminum alloy flip guardrail is convenient for patients to get on and off the stretcher, and is equipped with two safety belts.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue