

SKB2C01

HELICOPTOR STRETCHER



SKB2C01

FEATURES AND HIGHLIGHTS



Technical parameters

| | |
|-------------------------|-----------------|
| ▪ External Size (LxWxH) | 2180×620×180 mm |
| ▪ Self Weight | 18kg |
| ▪ Load Bearing | 300kg |



Technical configuration

| | |
|--|------|
| <input checked="" type="checkbox"/> Fix Belt | 3pcs |
| <input checked="" type="checkbox"/> Hanging Belt | 1pc |
| <input checked="" type="checkbox"/> Foot Fixator | 1pc |

Material

The body is made of non-toxic and pollution-free ABS material, which is light, pressure-resistant, easy to clean, fire-resistant, wear-resistant and corrosion-resistant.



Design

The structure of the stretcher focuses on the universality, flexibility and particularity of first aid, such as rescue in rugged mountains, air or sea. The frame is strong and durable, and the simple and reliable device enables the operator to take emergency measures safely and quickly.



Handle

The structure is made of reinforced aluminum alloy tube, with many arc-shaped smooth handles, which is convenient for many people to lift up in the process of use; metal fasteners are set on both sides of the stretcher, which can be used together with nylon fasteners.



Function

The suspension hook can be connected with the hook on the aircraft. Specially designed for field rescue.



Foot Fixator

The stretcher is equipped with adjustable foot fixator, which can be adjusted according to the height of the wounded.



Fix Belt

This product is equipped with three bandages as standard.



Hanging Belt

The optional kit includes neck bracket, head immobilizer.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue