



SKE087

RECLINING CHAIR



SAIKANG/Medical

Healthier
your life

SKE087

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Seat Size (LxWxH)	820×850×1070 mm
▪ Bed Size (LxWxH)	820×1600×450 mm


Backrest

Sectional type backrest gives the patient strong support and effectively protects the patient's spine.



Armrest

Solid wood armrest with rounded edges and a comfortable feel.



Cup Base

It is convenient for the patient to put the cup.



Design

Through the button to achieve a key flat, seconds into a recliner, can lie down to sleep.



Leather

Super strength PU material, Environmental protection.



Sponge

Flame retardant, 4 cm thickness, high-density foam, density 30 kg/m³.



Rotary Table

Solid wood rotary table, smooth and comfortable surface.



Castors

High quality medical silent casters, 3" casters, anti-static electricity, protect the ground, waterproof and anti-bacteria.




ULTRAFIBER LEATHER COLOR TABLE



Off white 5509-1



Off white 5509-2



Beige 5509-3



Pink 5510-45



Light rose red 2892-15



Rose red 2201-30



Wine red 2209-22



Light blue 5510-43



Sky blue 5510-7



Dark lake blue 2201-40



Light grey 5509-14



Gray 5510-12



Lake blue 5509-27



Navy blue 2201-42



Light green 2210-88



Green 5510-20

Note: The above is the physical scan , there may be a certain color difference in kind , are to prevail in kind.



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 · 24/7
 sales@chinasaikang.com
 www.saikangmedical.com



Company video



General catalogue