



SKE090

BLOOD DONATION CHAIR



SAIKANG/Medical

Healthier
your life

SKE090

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Seat Size (LxWxH)	650×800×450 mm
---------------------	----------------

Technical configuration

☒ Arm	2pcs
☒ Anti-slip Feet Cover	4pcs



Design

Scientific design, ergonomic, simple structure, convenient operation, applied to blood collection and rest.



Chair Surface

The PU leather surface is filled with high-quality sponge, which gives blood donors a sense of security at all times and a sense of comfort in the state of blood donation.



Material

Plastic steel frame, strong stability, durable use, beautiful appearance, comfortable texture, convenient cleaning.



Armrest

The armrest is filled with high-quality sponge, which can be rotated 360 degrees, which is convenient for blood donors to adjust, safe and reliable, and easy to operate.



Cabinet

There is a small cabinet next to the blood donors.



Foot Pad

Beautiful appearance, easy to use and has the functions of anti-slip, anti-collision, anti-shock, temperature resistance, noise reduction, etc.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue