



# SKH018

## ANESTHESIA TROLLEY

Registration code: SKM-B



**SAIKANG** Medical

Healthier  
your life

# SKH018

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 600×400×850 mm

### Technical configuration

- ☒ 4" Covered Castors 4pcs



#### Material

Made of the high-strength stainless steel, resistance of water, acid, alkali, fire and easy clean, not easily polluted by reagents, blood, drugs etc., satisfy the daily disinfectant and infection control with long life.



#### Table Top

Guardrail on three side to prevent the falling during moving.



#### Drawer

Upper drawer with plastic handle prevents cutting, one lock is equipped with two keys.



#### Cabinet

Lower cabinet with large capacity, plastic handle with lock, easy to storage large volume items.



#### Castors

4 nylon wheels Ø 100 mm with non-marking tires and plain bearings, two with brakes.





## **CONSTANT INNOVATION**



Design



Research



Manufacturing



Maintenance



Service

## **WORLDWIDE BRANCHES**

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

**SAIKANG** Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
 +86 18021231901 · 24/7  
 sales@chinasaikang.com  
 www.saikangmedical.com



Company video



General catalogue