



# SKH046-14

OVERBED TABLE




**SAIKANG** / Medical

Healthier  
your life

# SKH046-14

## FEATURES AND HIGHLIGHTS



### Technical parameters

External size (LxWxH) 890×400×895-1195mm

### Technical configuration

- |  |     |
|--|-----|
| <input checked="" type="checkbox"/> Air Spring for Height Adjustment | 1pc |
| <input checked="" type="checkbox"/> Phenolic Resin Table Top         | 1pc |

#### Material

U type base is made of powder coated steel, can be used in combination with bedside cabinets



#### Table top

Table top is made of HDPE, easy clean, high quality environmental material, light weight, corrosion resistance, impact resistance, easy to clean and disinfect, durable for use. With groove to put cups and pen.



#### Height adjustment

Height adjustment by gas spring, used single-handed according to the convenience of the user.



#### Casters

Easy to install, more stable and firmer. High quality medical caster, 1.5 inch / 2 inch mute wheel.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue