



SKR-AT625-1

ANESTHESIA TROLLEY

Registration code: SKM-B



SAIKANG Medical

Healthier
your life

SKR-AT625-1

FEATURES AND HIGHLIGHTS

Technical parameters


▪ Size (LxWxH)	625×475×930 mm
----------------	----------------

Technical configuration

☒ Sliding Side Shelf	1pc
☒ Dust Baskets	1pc
☒ Utility Container	1pc
☒ Central Lock	1set
☒ Trash Can	2pcs
☒ Anesthesia Frame	1pc

Drawers and Cabinet

Has two drawers, different size.
With doors, can put more things.



File Box

Side attached file box, for storage of
paper work.



Drawer Inside Dividers

With abs inserts inside, which can be
separated freely as required.



Drawer Card Holder

card insert from downside, easy keep
clean.



Anesthesia Frame

Anesthesia frame inserted into main
body column.



Powder Coated Column

Make the structure of the trolley more
stable and firm.



Table Top

Stainless steel side rail, avoid possible
fall down during movement of trolley.



Sliding Tray

Laterally extendable worktop.



Debris Basket

It is easy to put some files about
patient records.



Trash Can

There are two ABS dirt buckets of
the car body, different colors can
distinguish different kinds of dirt.



Castor

4 nylon wheels Ø100 mm with non-
marking tires and plain bearings, two
with brakes.



Heavy Duty Construction

Wheels assembled on the same
vertical line as trolley four columns,
more stable.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue