



SKR-CT232

CLINICAL TROLLEY

Registration code: SKM-B



SAIKANG Medical

Healthier
your life

SKR-CT232

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Size (LxWxH)	780×475×1000 mm
----------------	-----------------

Technical configuration

☒ File Box	1pc
☒ Plastic Basket	2pcs
☒ Sharps Boxes	2pcs
☒ Trash Can	2pcs
☒ Castor	4pcs



Material

Steel and aluminum structure inside, ABS materials outside, easy to clean.



Table Top

Guardrails on top of trolley is made of stainless steel, table top with soft glass, avoid possible movement and fall down during movement of trolley.



Push Handle

With push handle on both side, convenient for pushing easily.



Drawer

1 layer small drawer, 2-3 layer middle drawer, can put different sizes of medicine boxes and other medical materials, and equipped with central control lock.



Storage Shelf

Storage box which is designed fixed on the frame, large volume, easy to store and take out items, easy to clean.



Drawer Inside Dividers

With abs inserts inside, which can be separated freely as required.



Drawer Card Holder

Card insert from downside, easy keep clean.



File Box

Side attached file box, for storage of paper work.



Plastic Basket

With two plastic basket, there are two 2-liter Sharps boxes in the basket.



Trash Can

Two trash cans in different color with lid, controlled by foot switch.



Castor

4 nylon wheels Ø100 mm with non-marking tires and plain bearings, two with brakes.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue