



SKR-CT625-1

CLINICAL TROLLEY

Registration code: SKM-B



SAIKANG Medical

Healthier
your life

SKR-CT625-1

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Size (LxWxH)	625×475×920 mm
----------------	----------------

Technical configuration

☒ Trash Can	2pcs
☒ Castor	4pcs



Drawers

With Two shelves, can put more things, with one drawers, can put the medicine.



Table Top

Stainless steel side rail, avoid possible fall down during movement of trolley.



Trash Can

There are two ABS dirt buckets of the car body, different colors can distinguish different kinds of dirt.



Heavy Duty Construction

Wheels assembled on the same vertical line as trolley four columns, more stable.



Castor

4 nylon wheels Ø100 mm with non-marking tires and plain bearings, two with brakes.



Drawer Inside Dividers

With abs inserts inside, which can be separated freely as required.



Drawer Card Holder

card insert from downside, easy keep clean.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue