



SKR-CT670

CLINICAL TROLLEY

Registration code: SKM-B



SAIKANG Medical

Healthier
your life

SKR-CT670

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 670×450×825-1000 mm

Technical configuration

☒ Top Guardrail	Stainless Steel
☒ Trash Can	2pcs
☒ Castor	4pcs



Material

Made of the high quality cold rolled steel.
Color and size could be customized.
Stainless steel trolley column.



Drawer

2 small drawers with aluminum alloy handle, with PVC soft glass in drawer.



Castor

4 inch luxury noiseless covered casters, TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts.



Table Top

Stainless steel side rail, avoid possible fall down during movement of trolley. With stainless steel handrails on one side, convenient for pushing easily.



Trash Can

2 plastic removable trash can.





CONSTANT
INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE
BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue