



SKR-IB00

SIMPLE LAPTOP CART

Registration code: SKM-A



SAIKANG/Medical

Healthier
your life

SKR-IB00

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 533×455×750-1150 mm

Technical configuration

☒ Big Basket	1pc
☒ Tray	1pc
☒ Castor	4pcs



Column

Aluminum alloy column, the top equipment mesa.



Table Top

ABS material top, easy to clean and durable to scratch. Cover is free to open, laptop can be put in easily.



Compartment

Compartment can storage any laptop size under 17", is underneath top cover and lockable with key.



Holder for Mouse,

Double sided telescopic holder for mouse, left or right hand operating is both available.



Base

Base is ABS covered high quality metal structure, featured with both stability and easy-maintenance.



Castor

Casters is using medical silent material, double sided casters, stable and noiseless. Front 2 casters with brake, lockable at 360 degree, rear 2 casters free and swivel to any direction.



Big Basket

With basket on the back, easy to place goods, Environmentally plastic is non-toxic and easy to clean.



Option

Option configuration of high performance li-battery system, with LED Keyboard light, constantly providing stable power source, safe and convenient.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue