



SKR-J21

ENDOSCOPIC SYSTEM EQUIPMENT CART

Registration code: SKM-A



SAIKANG Medical

Healthier
your life

SKR-J21

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 665×630×1650 mm


Technical configuration

- ☒ Castor 4pcs



Column

Metal structure material comes from Baosteel from Fortune Global 500.



Holder of Screen

Above the table is equipped with a movable universal arm, which can hang the display, can rotate 360 degrees and adjust the height.



Working Table

Countertop height: 1125mm, dimensions: 527 (L) mm × 420mm (W), with a drawer under the countertop.



Movable Shelf

The height and number of each floor can be adjusted freely. The standard product is 2 adjustable movable shelves, which can be optionally added. The recommended load of the movable shelves is ≤ 20 kg.



Bumper Cover

With bumper, protect the cart during transportation.



Castor

4pcs of 100mm double side casters, each with brake.





CONSTANT
INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE
BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue