



SKR021-8

CASE HISTORY TROLLEY

Registration code: SKM-B



SAIKANG Medical

Healthier
your life

SKR021-8

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 750×410×950 mm

Technical configuration

- ☒ Drawer 2pcs
- ☒ Castor 4pcs

Material

Body material made of ABS; Powder coated steel column bearing.



Table Top

ABS injection moulded top, with grib and higher sided to prevent items flipping. With armrests on both sides.



Central Lock

Storage printed with label for catagory. With a central lock. With 50 holders for record cases.



Castor

Luxrury silent casters, with 2 lockable casters, stoppable at anytime, robust structure, universal direction casters.



Drawer

With 1 drawer, out height 115mm, inner dimension is 213×312×85mm.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 · 24/7
 sales@chinasaikang.com
 www.saikangmedical.com



Company video



General catalogue