



X01-5

BABY CRIB

Registration code: SK-B3

X01-5

FEATURES AND HIGHLIGHTS



Technical parameters

- External size (LxWxH) 810×520×990 mm
- Trendelenburg 18°

Technical configuration

- Bassinet 1pc
- 3" Covered Castors 4pcs
- IV Pole Prevision 1pc
- IV Pole 4pcs
- Utility Basket 1pc
- Adjustable Air Spring 4pcs



Material

Stainless steel frame, baosteel from Fortune Global 500.



Painting

11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Baby Basin

Transparent baby basin, environmentally acylic non-toxic, odorless, can be sterilized at temperatures above 100 °C.



Trendelenburg

The backward tilt design adopts pneumatic lifting mode, which is stable and can prevent the infant from vomiting after eating.



Castor

3" covered casters. TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts.



Utility Basket

There is a litter basket at the bottom to store the baby's items.




Mattress

Custom mattress, easy disassembly and cleaning, optional.



IV Pole

Full SS IV pole, bearing 15kg.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue