



X01

BABY CRIB

Registration code: SK-B3



SAIKANG / Medical

Healthier
your life

X01

FEATURES AND HIGHLIGHTS




Technical parameters

- External size (LxWxH) 810×520×990 mm
- Trendelenburg 18°

Technical configuration

- Transparent baby basin 1pc
- 3" Castors 4pcs
- IV pole 1pc
- Dustproof 4pcs



Material

Stainless steel frame, baosteel from Fortune Global 500.



Mattress

Custom mattress, soft and comfortable material, dustproof, antibacterial and easy to clean.



Trendelenburg

Controlled by mechanical structure to achieve trendelenburg, convenient and safe.



Baby Basin

Transparent baby basin, environmentally acyclic non-toxic, odorless, can be sterilized at temperatures above 100 °C.



Castor

3" Bulging pole wheel, two with brakes, easy to install, stable and smooth.




IV Pole

Full SS IV pole, bearing 15kg.



Dustproof

4 dust plugs, effectively prevent dust, detachable and easy to clean.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue