



X10

EXAMINATION TABLE


Registration code: SK-B1

SAIKANG / Medical

Healthier
your life

X10

FEATURES AND HIGHLIGHTS



Technical parameters

- | | |
|-------------------------|-----------------|
| ▪ External size (LxWxH) | 1900x600x680 mm |
| ▪ Back-rest adjustment | 0-65° (±10°) |


Technical configuration

- | | |
|--|-------|
| <input checked="" type="checkbox"/> Anti-slip bed feet cover | 4 pcs |
| <input checked="" type="checkbox"/> Mechanical level to adjust back rest | 1 set |



Leather

Yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance.



Sponge

6cm thickness, high-density foam, density 30 kg/m³.



Painting

11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Metal material

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth wedding.



Frame

1.2-1.5 mm thickness stainless steel frame, 250kg bearing loading.



Feet cover

4 pcs anti-slip bed feet covers





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue