



X20

PEDIATRIC EXAM TABLE

Registration code: SKD-A



SAIKANG / Medical

Healthier
your life

X20

FEATURES AND HIGHLIGHTS

Technical parameters

- External size (LxWxH) 1800x650x830 mm

Design

Double-sided cartoon painting, unique and novel.




Material

Multi-layer solid wood board thickness 18mm, surface painted (glue blue).



Door Locker

Double open door locker, size: 520×500mm.




Mattress

The mattress can be removed, the mattress size: 1100×610×60mm.



Leg Plate

The leg plate is under the bed board and can be pulled out 260mm. The leg pad size is 520×260×90mm.



Stair Pedal

Stair pedal design for children to get on and get out of bed.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue