



# CT2k

## CHILDREN BED

Registration code: SK-C2




**SAIKANG** / Medical

Healthier  
your life

# CT2k

## FEATURES AND HIGHLIGHTS




### Technical parameters

▪ External size (LxWxH)	1800x730x500 mm
▪ Mattress platform	1580x660 mm
▪ Safe working load	250 kg
▪ Back-rest adjustment	0-75° (±10°)
▪ Knee-rest adjustment	0-35° (±10°)

### Technical configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs
<input checked="" type="checkbox"/> Aluminium Alloy Side Rail with Extra Protection	1set
<input checked="" type="checkbox"/> Aluminium Alloy Bed Ends	1set
<input checked="" type="checkbox"/> Anti-slip Bed Feet Cover	4pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Shoes Holder	1pc



**Steels | Painting**

Baosteel from Fortune Global 500  
11 Process epoxy painting, ASTM  
testing anti-bacterial, paint thickness  
0.12mm, brightness 60°, paint can  
resist 50kg impact.



**Welding**

Panasonic Robotic ensure,  
360° full smooth wedding.



**Hand Control**

Retractable hand crank, 100,000  
times remaining intact.  
Extreme protection device, warning  
when tipped to the max position.



**Head & Foot Board**

Aluminum alloy flat tube, MDF  
material fire-proof plate. size of the  
head&foot board can be customized.  
With carton photograph.




**Siderail**

Aluminum alloy foldable guardrail  
with protection handle, modern  
and beautiful, easy to install, easy to  
clean, anti-pinch, protect patient's  
safety.



**Fix**

4 pcs Anti-slip bed feet covers.



**Platform**

Punching platform, smooth surface  
no welding scar, bearing 250KG  
load. Option resin X-ray translucent  
platform.



**Drainage Hooks**



**Shoes holder**



**IV Pole**

Stainless Steel Column with plastic  
hook.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 📞 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue