



# SKR-ET850

## EMERGENCY TROLLEY



**SAIKANG** / Medical

Healthier  
your life

# SKR-ET850

## FEATURES AND HIGHLIGHTS

### Technical parameters

- External size (LxWxH) 850x520x990 mm
- Castors 125 mm




#### Material

It is mainly composed of aluminum, stainless steel and ABS engineering plastic structure.



#### Drawer handle

The color of the drawer handle can be changed, easy to assemble.



#### Defibrillator holder

Two-layer rotating defibrillator holder, can easily view and control the defibrillator at any Angle.



#### Drawer

Drawer guide rail adopt damping buffer rail.



#### Castor

4 inch luxurious noiseless castors, and two with brake high strength polyurethane material.



#### Drawer label card

Card insert from downside, easy keep clean. Color is also optional.



#### ABS inserts

With ABS inserts inside, which can be separated freely as required.



#### Laterally extendable worktop



#### Safe lock

Adopt C-level safe lock, make medicine storage safer.



#### Dust baskets

Two dust baskets are divided into two colors, convenient for medical staff to classify different types of garbage





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

+86 18021231901 · 24/7

sales@chinasaikang.com

www.saikangmedical.com



Company video



General catalogue