



# SKE001-12

SOFA CHAIR



**SAIKANG** / Medical

Healthier  
your life

# SKE001-12

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ External size	2100×850×620mm
▪ Seat size	830×850×420mm (backrest height: 960mm)
▪ Safe working load	200kg
▪ Color	MOQ>20 colors can be customized




#### Design

Ergonomic design to reduce muscle tension and stress, stylish appearance, leather armrests, high-resilience sponge inside, comfortable sitting and lying position, durable and easy to clean.




#### Frame

The internal frame is made of 1.0 thickness iron sheet folded square tube, safe load up to 200kg.



#### Surface material

PU leather, moisture and water proof, thickened sponge inside, comfortable and high wear resistance.



#### Bed legs

Bed legs are equipped with plastic cover to prevent possible noise caused by friction with the floor surface; when folded the backrest is securely fixed, to enhance patient's safety.




#### Caster

Equipped with 3-inch silent caster\*4pcs, 1.5-inch silent caster\*4pcs, smooth, easy operable.



#### Function

Telescopic springs inside, material 65Mn, self-positioning design, two easy switchable modes.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### South America Branch

**Saikang Medical Technology, Inc**  
W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

### East Africa Branch

**Saikang Medical Technology Ltd**  
Mezzanine Floor, Taj Tower,  
upper hill, nairobi, Kenya  
TEL: +254714084806

### South East Asia Branch

**Saikang Medical Technology Ltd**  
Sala Thammassop Subdistrict  
Thawi  
Watthana District  
Bangkok Thailand  
TEL: +66 642741046

### UK Branch

**Uk Healthward International Limited**  
DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +31 (0) 202111106

### Russia Branch

**LLC 'SAIKANG'**  
454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya Naberezhnaya,  
Russia  
TEL: +7 (919)300-00-06

### Mid East Branch

**Saikang Medical Supplies (Mid East) LLC**  
Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ export@saikangmedical.com

🌐 www.saikangmedical.com



Company video



General catalogue