



SKE001-13

ACCOMPANY CHAIR



SAIKANG / Medical

Healthier
your life

SKE001-13

FEATURES AND HIGHLIGHTS

Technical parameters

- External size 650×2060×450mm
- Seat size 780×650×450mm(backrest height: 1060mm)
- Safe working load 200kg
- Color MOQ>20 colors can be customized




Design

The overall design is simple yet high-grade, fold into a chair mode during the day and unfold into a bed mode at night, safe and reliable. Elegant curved armrests, ergonomic design, stylish and comfortable.




Frame

The internal frame is made of 1.0 thickness iron sheet folded square tube, safe load up to 200kg.




Surface material

PU leather, moisture and water proof, thickened sponge inside, comfortable and high wear resistance.



Armrest

Solid wood armrests provide a luxurious and gentle touch.



Caster

Equipped with 3-inch silent caster*4pcs, 1.5-inch silent caster*4pcs, smooth, easy operable.



Function

Telescopic springs inside, material 65Mn, self-positioning design, two easy switchable modes.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc
W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd
Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd
Sala Thammassop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

UK Branch

Uk Healthward International Limited
DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG'
454004, Russian Federation,
Chelyabinsk,
Universitetskaya Naberezhnaya,
Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC
Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ export@saikangmedical.com

🌐 www.saikangmedical.com



Company video



General catalogue