




# CQ8k

## ELECTRIC CHILDREN BED

Registration code: SKD-C





**SAIKANG** / Medical

Healthier  
your life

# CQ8k

## FEATURES AND HIGHLIGHTS



### Technical parameters

External size (LxWxH):	1975x880x465-760 mm
Mattress platform:	1730x780 mm
Safe working load:	200 kg
Back-rest adjustment:	0-75° (±5°)
Knee-rest adjustment:	0-35° (±5°)
Tredelenburg:	0-12°
Reverse Trendelenburg:	0-12°

### Technical configuration

<input checked="" type="checkbox"/> Linear motor	4pcs
<input checked="" type="checkbox"/> Stainless steel side rail with extra protection	1set
<input checked="" type="checkbox"/> PP bed ends with safe lock	1set
<input checked="" type="checkbox"/> PP bed platform, removeable, washable	5pcs
<input checked="" type="checkbox"/> 5" double side cators with central locking system	4pcs
<input checked="" type="checkbox"/> IV pole	1set
<input checked="" type="checkbox"/> IV pole prevision	4pcs
<input checked="" type="checkbox"/> Drainage hook	2pcs

### Optional

- Trapeze pole
- Directional wheels function

# QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



## Design

This Bed uses unique design to meet hospital strict gremfree request. Head & Foot board and platform are all seamless, which supports effective sterilization.



## Steels | welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



## Plastic | Painting

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa.  
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



## Electric Motor

American motor, with UL/ROHS/EN standard. Noise under 50dB, with lifetime 20,000 times movement.

## Accessories



## Platform

Five PP material bed platforms are washable and easily movable.



## Castors | Brake system

Castors with TPR tire, no worn-out after running 30KM, united forming without bolts.  
Pass Dynamic test: bearing 120kg run 30KM and pass obstacles 500time.  
Brake uses 301 stainless steel. Castors' brakes keep qualified after 10,000 times test.

**Head & Foot Board**

Modern design PP head&foot board.  
Easy to be cleaned and installed.



**Handle Controller**

Handle controller contains flexible spiral cable for easy mobility and operation in optimal position. Patient and nurse can control all necessary functions.



**Siderail**

Side rail made from stainless steel and take 12-staff design. And it with stainless steel protection for patients can get on or off bed easily.




**Drainage Hooks**

Mobile plastic drainage hooks.



**IV Pole**

Full SS IV pole, bearing 2.5kg per hook.




Safe working load

**200**kg




Back-rest adjustment

**75°**




Knee-rest adjustment

**35°**




Height adjustment

**465-760**mm



Reverse trendelenburg

**12°**



Trendelenburg

**12°**





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### South America Branch

**Saikang Medical Technology, Inc**  
 W3400 Corrientes, Argentina  
 TEL: +5493794312231  
 Av. Senador Vergueiro, 4204 - sala 31  
 Rudge Ramos - São Bernardo do Campo - SP, Brazil  
 TEL: +55 (11) 2758-0438

### East Africa Branch

**Saikang Medical Technology Ltd**  
 Mezzanine Floor, Taj Tower,  
 upper hill, nairobi, Kenya  
 TEL: +254714084806

### South East Asia Branch

**Saikang Medical Technology Ltd**  
 Sala Thammassop Subdistrict  
 Thawi  
 Watthana District  
 Bangkok Thailand  
 TEL: +66 642741046

### UK Branch

**Uk Healthward International Limited**  
 DEPT 111,  
 196 High Road Wood Green,  
 London, N22 8HH  
 United Kingdom  
 TEL: +31 (0) 202111106

### Russia Branch

**LLC 'SAIKANG'**  
 454004, Russian Federation,  
 Chelyabinsk,  
 Universitetskaya Naberezhnaya,  
 Russia  
 TEL: +7 (919)300-00-06

### Mid East Branch

**Saikang Medical Supplies (Mid East) LLC**  
 Musaffah  
 P.O Box 53200, Abu Dhabi,  
 UAE  
 Tel: +971 585710652

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
 STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ export@saikangmedical.com

🌐 www.saikangmedical.com



Company video



General catalogue