



# A100-4A

## ELECTRIC OPERATION TABLE

Registration code: SKL-A



**SAIKANG** / Medical


Healthier  
your life

# A100-4A

## FEATURES AND HIGHLIGHTS

### Technical parameters

External Size (LxWxH):	2060x500x710-1010 mm
Built-in Kidney Bridge:	0-110 mm
Back Plate Up/Down:	80°/20°
Head Plate Up/Down:	40°/90°
Leg Plate Up/Down:	15°/90°
Trendelenburg/Reverse Trendelenburg:	25°/25°
Lateral Tilt Left/Right:	15°/15°
Horizontal Sliding:	350 mm
Battery Standby Time:	80 h
Power Supply:	220V 50HZ 350VA
Power Control Circuit:	12V
Plug:	as required
X-rays translucent:	YES



- C-arm compatible: YES
- Multi-functional orthopedic frame: YES

### Technical configuration

<input checked="" type="checkbox"/> Bed Frame	1pc
<input checked="" type="checkbox"/> Anesthesia Screen Frame	1pc
<input checked="" type="checkbox"/> Mattress	1set
<input checked="" type="checkbox"/> Kidney Bridge Handle	1pc
<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Shoulder Support with Cushion	2pcs
<input checked="" type="checkbox"/> Body Support with Cushion	2pcs
<input checked="" type="checkbox"/> Foot Support with Cushion	2pcs
<input checked="" type="checkbox"/> Hand Controller	1pc

#### Remote Panel & Column Panel

Leg Plate With Imported Gas Spring For Easy Control.



#### Kidney Bridge

Built-in kidney bridge, more convenient and nice appearance.



#### Wheel

6" castors import from germany.



#### Leg Plate

Leg Plate With Imported Gas Spring For Easy Control.



#### Material

The base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.



#### Braking System

Pedal Brake, Easy To Move And Lock.



#### Powder Supply

Internal battery maintaining power supply for 50-80 operations and with low battery warning function.



#### Shoulder Support with Cushion



Body Support with Cushion



Arm Support with Cushion



Leg Support with Cushion



Anesthesia Screen Frame



Lowest Position of Table Top  
**710mm**



Highest Position of Table Top  
**1010mm**



Back Plate Up  
**0°-80°**



Back Plate Down  
**0°-20°**



Reversed Trendelenburg  
**0°-25°**



Trendelenburg  
**0°-25°**



Lateral Tilt Left  
**0°-15°**



Lateral Tilt Right  
**0°-15°**



Horizontal Sliding  
**350mm**



Head Plate Up  
**0°-40°**



Head Plate Down  
**0°-90°**



Leg Plate Rest  
**-90°~15°**





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### South America Branch

#### Saikang Medical Technology, Inc

W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

### East Africa Branch

#### Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,  
upper hill, nairobi, Kenya  
TEL: +254714084806

### South East Asia Branch

#### Saikang Medical Technology Ltd

Sala Thammasop Subdistrict  
Thawi  
Watthana District  
Bangkok Thailand  
TEL: +66 642741046

### UK Branch

#### Uk Healthward International Limited

DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +31 (0) 202111106

### Russia Branch

#### LLC 'SAIKANG'

454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya  
Naberezhnaya, Russia  
TEL: +7 (919)300-00-06

### Mid East Branch

#### Saikang Medical Supplies (Mid East) LLC

Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@saikangmedical.com

🌐 www.saikangmedical.com



Company video



General catalogue