



# Q2w8c

## MANUAL DIALYSIS BED



**SAIKANG** Medical


Healthier  
your life

## Q2w8c

### FEATURES AND HIGHLIGHTS

#### Technical parameters

▪ External size (LxWxH)	2185x880x540mm±10 mm
▪ Mattress platform	1940x780mm±10mm
▪ Back-rest adjustment	0-75° (±10°)
▪ Knee-rest adjustment	0-35° (±10°)




#### Technical configuration

☑ Manual Cranks	2pcs
☑ Aluminium Alloy Side Rail with Extra Protection	1 set
☑ 5" Double side castors, Central Locking System	4pcs
☑ PP bed ends with safe lock	1 set
☑ PP bed platform, Removeable & Washable	4 pcs
☑ Iv pole	1pc
☑ Iv pole prevision	4pcs



#### Design

Modern design. Use new and strong inserted PP fire-proof panel avoid off. Easy to be cleaned and installed



#### Mechanical functions

2 Functions: Backrest adjustment, Legrest adjustment



#### Castors

5" Double side castors with central locking system  
Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500 times



#### Platform

PP platform, removeable, washable




**Head & Foot Board**

Modern design. Use new and strong inserted ABS fire-proof panel avoid off. Easy to be cleaned and installed.

**Siderail**

L-shaped foldable guardrail (Patent), bended tube, solid installing, 0 shaking, anti-pinch, protect ensure patient's safety

**Hand Cranks**

The bed with two functions, controlled by manual cranks. Retractable hand crank, 100,000 times remaining intact

**Mattress Stopper**

Platform integrated mattress stopper on both sides of the back and leg section.

**Utility Shelf**

Detachable utility shelf for patients' personal stuff.

**IV Pole**

Full SS IV pole, bearing 2.5kg per hook





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### South America Branch

**Saikang Medical Technology, Inc**  
W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

### East Africa Branch

**Saikang Medical Technology Ltd**  
Mezzanine Floor, Taj Tower,  
upper hill, nairobi, Kenya  
TEL: +254714084806

### South East Asia Branch

**Saikang Medical Technology Ltd**  
Sala Thammasop Subdistrict  
Thawi  
Watthana District  
Bangkok Thailand  
TEL: +66 642741046

### UK Branch

**Uk Healthward International Limited**  
DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +31 (0) 202111106

### Russia Branch

**LLC 'SAIKANG'**  
454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya  
Naberezhnaya, Russia  
TEL: +7 (919)300-00-06

### Mid East Branch

**Saikang Medical Supplies (Mid East) LLC**  
Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

# SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎ +86 18021231901 · 🕒 24/7  
✉ export@saikangmedical.com  
🌐 www.saikangmedical.com



Company video



General catalogue