



A100-4

ELECTRIC OPERATING TABLE


Registration code: SKL-A

SAIKANG / Medical

Healthier
your life

A100-4

FEATURES AND HIGHLIGHTS



Technical parameters

External Size (LxWxH):	2070x550x700-1000 mm
Built-in Kidney Bridge:	0-120 mm
Back Plate Up/Down:	75°/20°
Head Plate Up/Down:	45°/90°
Leg Plate Up/Down:	15°/90°
Trendelenburg/Reverse Trendelenburg:	25°/25°
Lateral Tilt Left/Right:	15°/15°
Horizontal Sliding:	300 mm
X-rays translucent:	YES
C-arm compatible:	YES
Multi-functional orthopedic frame:	YES

Technical configuration

<input checked="" type="checkbox"/> Bed Frame	1pc
<input checked="" type="checkbox"/> Anesthesia Screen Frame	1pc
<input checked="" type="checkbox"/> Mattress	1set
<input checked="" type="checkbox"/> Kidney Bridge Handle	1pc
<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Shoulder Support with Cushion	2pcs
<input checked="" type="checkbox"/> Body Support with Cushion	2pcs
<input checked="" type="checkbox"/> Foot Support with Cushion	2pcs

Remote Panel

All function controlled by hand controller, and can be locked.



C-arm

Extra-long horizontal sliding, available for c-arm operation.




Kidney Bridge

Built-in kidney bridge, more convenient and nice appearance.



Wash Basin Hole

Five section table top including Head, Back, Shoulder, Hip and Leg, all section applicable for X ray use and under each section with X ray cassette layer.



Head Plate

Detachable double head plate and leg plate meet operation position needs and use Germany cylinder assist system.



Powder Supply


Internal battery maintaining power supply for 50-80 operations and with low battery warning function.



Shoulder Support with Cushion



Body Support with Cushion



Arm Support with Cushion



Leg Support with Cushion



Lowest Position of Table Top
800mm



Highest Position of Table Top
1000mm



Longitudinal Adjustment
0-300mm



Longitudinal Adjustment
0-300mm



Reversed Trendelenburg
0°-22°



Trendelenburg
0°-25°



Lateral Tilt Left
0°-20°



Lateral Tilt Right
0°-20°



Kidney Bridge Elevation
110mm



Anorectum Operation
Position



E.N.T Operation Position



(Optional) Head Rest



(Optional) Lateral Arm
Support



(Optional) Waste Basin For
Gynecology



(Optional) Orthopedics Traction
Frame





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina
 TEL: +55 (11) 2758-0438
 Rua Universal, Hollywood,
 São Bernardo do Campo,
 Brazil
 TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,
 upper hill, nairobi, Kenya
 TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict
 Thawi
 Watthana District
 Bangkok Thailand
 TEL: +66 642741046

UK Branch

Uk Healthward International Limited

DEPT 111,
 196 High Road Wood Green,
 London, N22 8HH
 United Kingdom
 TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG'

454004, Russian Federation,
 Chelyabinsk,
 Universitetskaya
 Naberezhnaya, Russia
 TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah
 P.O Box 53200, Abu Dhabi,
 UAE
 Tel: +971 585710652

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
 STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@saikangmedical.com

🌐 www.saikangmedical.com



Company video



General catalogue