



A98-4

ELECTRIC OBSTETRIC BED




SAIKANG Medical

Healthier
your life

A98-4

FEATURES AND HIGHLIGHTS




Technical parameters

External Size (LxW)	1980x960mm(±30mm)
Height	650-940mm(±20mm)
Back-rest Adjustment	0-60°
Reverse Trendelenburg	0-7°
Leg-rest Adjustment	-9°-35°, leg section is removable
Weight capacity	200KG



Technical configuration

<input checked="" type="checkbox"/> Linear motor	3pcs
<input checked="" type="checkbox"/> 5" Double side castors, central locking	4pcs
<input checked="" type="checkbox"/> Central locking pedal	2pcs
<input checked="" type="checkbox"/> ABS head board with safe lock	1pc
<input checked="" type="checkbox"/> ABS side rail with angle indicator	1set
<input checked="" type="checkbox"/> Hand grip	1pair
<input checked="" type="checkbox"/> Leg holder	1pair
<input checked="" type="checkbox"/> Basin	1pc
<input checked="" type="checkbox"/> Hand controller	1pc
<input checked="" type="checkbox"/> Power line	1pc
<input checked="" type="checkbox"/> Leg strip	1 pair
<input checked="" type="checkbox"/> IV pole prevision	2pcs



Material

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.

**Bed surface**

Bed surface by plastic steel, made by virgin material. The cover is detachable, and easy to clean.

**Head Board and Guardrail**


The Head board and side rail made by virgin ABS material. The head board is movable.

**Hand controller**

Equipped with one hand controller.

**Handle**

The hand bar can be dismantled and movable.

**Lithotomy rod**

The leg support can be dismantled and movable.

**Caster**

5 inch central controlled caster for moving obstetric table. TPR tire, quiet and wear-resistant, can be washed and sterilized by high temperature, save anti-winding hard shell, united forming without bolts.

**Accessory**

Equipped with stainless steel sewage basin and other various accessories.

**Mattress**

The mattress adopts high density memory sponge and is once formed. Mattress is movable.

**Foot section**

Foot section can be retracted, convenient for gynecological and delivery surgical.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc
 W3400 Corrientes, Argentina
 TEL: +55 (11) 2758-0438
 Rua Universal, Hollywood,
 São Bernardo do Campo,
 Brazil
 TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd
 Mezzanine Floor, Taj Tower,
 upper hill, nairobi, Kenya
 TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd
 Sala Thammassop Subdistrict
 Thawi
 Watthana District
 Bangkok Thailand
 TEL: +66 642741046

UK Branch

Uk Healthward International Limited
 DEPT 111,
 196 High Road Wood Green,
 London, N22 8HH
 United Kingdom
 TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG'
 454004, Russian Federation,
 Chelyabinsk,
 Universitetskaya
 Naberezhnaya, Russia
 TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC
 Musaffah
 P.O Box 53200, Abu Dhabi,
 UAE
 Tel: +971 585710652

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
 STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 📞 +86 18021231901 · 🕒 24/7
 ✉️ export@saikangmedical.com
 🌐 www.saikangmedical.com



Company video



General catalogue