



SKE910

PEDIATRIC INFUSION CHAIR



SKE910


FEATURES AND HIGHLIGHTS

Design: Cartoon-style design, tailored to childrens' psychological needs. This can help children relax and better receive the treatment. The independent seat design helps to avoid overcrowding and prevent cross-infection.

Material: Use European E0 grade standard particleboard, made by a top domestic manufacturer "Furen". It has a static bending strength of 13.7MPa, moisture content of 6.8%, density of 0.67g/cm³, and formaldehyde emission of 0.047mg/m³. It is durable, sturdy, and environmentally friendly.

Technical parameters

Overall Dimensions	1236×900×810cm
Seat Size (LxW)	1180×350mm



Veneer

The board is covered with "Xiate" decorative paper, which is wear-resistant and not easy to change color. The surface is smooth and glossy, with a good overall effect.



Edge

Use high-quality solid wood particle board as the base and PVC as edges. "Lisheng" hot melt glue is used to press the edges at a high temperature of 200°C to prevent moisture from penetrating. The beveled edge is sealed with high-quality kraft paper.



Upholstery

The fabric has a dry rubbing resistance level of 4 or higher, anti-static and anti-slip. The seat cushion is made of high-resilience sponge, which is comfortable for long-time sitting.



Chair leg

L-shaped chair leg design, stable and strong load-bearing.




Seat

The drawer-style stool has storage function, which saves space for parents to accompany better.



Side board

Heart-shaped side plate adds a touch of playfulness, helping pediatric patients relax for easier treatment.



Corner

Rounded corners to prevent collision and accidental injury.



Seat

Large seating area provides enough space for family accompany and toys.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc
 W3400 Corrientes, Argentina
 TEL: +55 (11) 2758-0438
 Rua Universal, Hollywood,
 São Bernardo do Campo,
 Brazil
 TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd
 Mezzanine Floor, Taj Tower,
 upper hill, nairobi, Kenya
 TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd
 Sala Thammassop Subdistrict
 Thawi
 Watthana District
 Bangkok Thailand
 TEL: +66 642741046

UK Branch

Uk Healthward International Limited
 DEPT 111,
 196 High Road Wood Green,
 London, N22 8HH
 United Kingdom
 TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG'
 454004, Russian Federation,
 Chelyabinsk,
 Universitetskaya
 Naberezhnaya, Russia
 TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC
 Musaffah
 P.O Box 53200, Abu Dhabi,
 UAE
 Tel: +971 585710652

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
 STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 📞 +86 18021231901 · 🕒 24/7
 ✉️ export@saikangmedical.com
 🌐 www.saikangmedical.com



Company video



General catalogue